

Product datasheet for **RN210434**

Tspan3 (NM_001005547) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tspan3 (NM_001005547) Rat Untagged Clone
Tag: Tag Free
Symbol: Tspan3
Synonyms: Tm4sf8
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN210434 representing NM_001005547
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGCCAATGCGGCATCACCTCCTCTAAGACCGTGTGGTCTTCTCAACCTCATCTTCTGGGGGGCGG
 CTGGCATCCTCTGCTACGTGGGAGCCTATGTCTTCATCACTTACGACGACTATGACCACTCTTTGAGGA
 TGTGTACACGCTCTCCCTGCCGTGGTATCATGGCTGTAGGGCCTTGCTTTCATCATCGGCCTGATC
 GGCTGCTGTGCTACAATCCGGGAGAGCCGCTGTGGGCTGGCCACGTTTGTCTTCATCCTGCTCTTGGTTT
 TCGTCACAGAAGTAGTAGTTGTGGTTTTGGGATATGTGTACAGAGCAAAGGTGGAGAATGAGGTTGATCG
 TAGCATTACAGAAAGTGATAAGACCTATAACGGGACCAACTCTGATGCTGCCAGCCGTGCTATTGACTAT
 GTGCAGAGACAGTTGCACTGCTGTGGAATTCATAACTATTCAGACTGGGAAAATACAGATTGGTTCAAAG
 AGACAAAAACCAGAGTGTCTCTTAGCTGCTGCAGAGAGACTGCCAGAAGCTGTAATGGCAGCCTGGC
 CAACCCCTCCGACCTCTATGCCGAGGGCTGCGAGGCTCTTGTGTGAAGAAGCTACAAGAAATCCTGATG
 CATGTTATCTGGCAGCTTTGGCCTTTGCAGCTATTCAGCTGCTGGGCATGCTGTGTGCATGCATCGTGC
 TGTGCAGACGGAGTAGAGACCCTGCCTATGAGCTCCTCATCACTGGTGGGACCTACGCATAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001005547
Insert Size: 762 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001005547.1</u> , <u>NP_001005547.1</u>
RefSeq Size:	1644 bp
RefSeq ORF:	762 bp
Locus ID:	300733
Cytogenetics:	8q24