

Product datasheet for **RN210450**

Ly6g6e (NM_001001972) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Ly6g6e (NM_001001972) Rat Untagged Clone

Tag: Tag Free

Symbol: Ly6g6e

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN210450 representing NM_001001972
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGGCTGTCCAGCGCCTTCCTTGGCCTCCTGTTCTCTCGGGGACACTGGGTCTCACCCTTCCCCAA
CCAGAGGACGACTCCACTGTTACACCTGCAGCTTTGCCAAGCCCTGCTACCCTGTTCTCACGGAGTGTCTG
GGAAGATGAGGTGTGCGGTGTCACTGTTGGCACCTCAGAGCAAAACAAGGATGTCATCGAGCGGAAAGGC
TGCTCTCCCAAGAGCCAGTGCCTCTGCTGGGCCATACTACCTACTGGTCACGGTCTACACTCTCCAGC
ATCAGTGCTGTGAGCAGGATCTGTGCAACACTGCTGCCTCACAGCGACCCCGAGCCTTCCCCTCATGAC
TCTGCTACCCCTTGACGCCATCTTTGGCTGGAGGTACATATCGTCCTC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001001972

Insert Size: 402 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001001972.1](#), [NP_001001972.1](#)

RefSeq Size: 959 bp

RefSeq ORF: 402 bp

Locus ID: 406866

Cytogenetics: 20p12

Gene Summary: may play a role in immune response [RGD, Feb 2006]