

Product datasheet for **RN210451**

Gpbar1 (NM_177936) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gpbar1 (NM_177936) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Gpbar1
Synonyms:	Tgr5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN210451 representing NM_177936 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATGTCACACAACACCACTGAGCTGTGAGCCATCCCAGAGGGGTTCCAGGAGCTTCCCTGGTCTGG
CAAGCCTCATCGTCATCGCCAACCTGCTCCTGGCCCTAGGCATTGTCTGGACCGCCACTTACGCAGCCC
ACCTGCTGGCTGCTTCTTTCTAAGCCTACTACTAGCTGGGCTACTCACAGGGTTGGCACTGCCACGCTG
CCTGGGCTATGGAATAGGAGCCATCAGGGTACTGGTCTGCCTCCTTCTCCACTGGCCCCCACTTTT
GTTTCTCTCCCTGCTTGCCAATCTGCTGCTGGTGCATGGGGAACGCTACATGGCAGTGTGCAGCCACT
CCGGCCCCATGGGAGTGTGCGGCTAGCCCTGTTCTCACCTGGATCAGCTCCCTGCTCTTGGCAGCCTG
CCTGCTCTGGGCTGGAACCACTGGAGTCTGGTGCCTGCACTGCCAGCTCCCAGGCTATCTTCCAGCCCCCT
ACCTTTACCTCGAAGTCTATGGGCTCCTGCTGCCCGCTGTGGGGCCACTGCCCTTCTCTGTCCGAGT
GTTGGCCACTGCCACCACCAGCTGCGGGAGATCCGCAGACTGGAGCGGGCGGTGTGCCGTGATGCACCC
TCAACCCTAGCGAGGGCTCTCACCTGGAGGCAGGCTAGGGCACAGGCAGGAGCCACACTGCTTTTTTGC
TGTGTTGGGGCCCTATGTGGCCACATTGCTCCTGTGAGTCTTGGCCTATGAGCGGGCCACCAGTGG
GCCTGTAACCTGTTATCTCTCATCTCATTGGGCAGTGCAGTGTGCAGTTGTGCTGTGGCCATGGGT
CTGGGTGATCAGCGCTACACGGCCCCCTGGAGGACAGCTGCCAAAGGTGGCTACAAGTGCTTCGAGGAA
GACCCAAGAGGGCCAATCCAGGCCCCAGCACTGCCTACCCTCCAGTAGCCAATGCAGCACTGACTTGA
CTTGAAT**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja1440_b04.zip



[View online >](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_177936
Insert Size:	990 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_177936.1</u> , <u>NP_808797.1</u>
RefSeq Size:	990 bp
RefSeq ORF:	990 bp
Locus ID:	338443
UniProt ID:	<u>Q80T02</u>
Cytogenetics:	9q33
Gene Summary:	a bile acid receptor [RGD, Feb 2006]