

Product datasheet for **RN210523**

Pdcd2l (NM_001109544) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pdcd2l (NM_001109544) Rat Untagged Clone
Tag: Tag Free
Symbol: Pdcd2l
Synonyms: RGD1308249
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN210523 representing NM_001109544
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGCCGTCAGGAAGCCGGTCTGCTGGGATTTTCGAGACTCGGCGGTGAAGGGCTGCCCAAGGGAC
CCGCCGCCTGGACGGCTAGCAAGATAGGGCGCCTCCGGACGCCCTGCCGGCCTTACTGCCCGGGCC
GCGGTGTGACCGCTGCGCACAGCCGCTCGCGCTCGTGCAGGTGTACTGCCCGCTGGAAGGCTCCCCG
TTCCATCGGCTTCTGTAGTTTTTCGCTGTGCCGTCGGGACTGCAGCGACAACCAGACGCGCAGCTGGA
AGGTGTTCCGCTCCCAGTGCCTGCAAGTGCCAGAGAAAGAGACTCGGAACGGACAGAATCAGAGTGAAGG
CCTTGCCGCTGAGAACTGGTGTGAGGGCTCACATGACTGGGAAGTGATACCGAGGAGACGCTCCGCCA
CCACCCAGCAGATTTGGGAAGGGATTCTAATGATGCCAGAGCGCTAGACTGGACTGACAAGCTCCAAG
CTCTTCACCTGCAGGATAGTGCCCTGGCTGTTGCCTGCCCTTCCCTTCAGTCCCACGGAAGTGCCACA
GTTCCAGCCCTATTATACAGTGTGCGAAGAGGATGAGTACGGCAACGTTGCCGATCTTGACCACGCC
CACAGCTCCTGCAGGAGTACCAGCAGAGAGAAGGTGTTGACATGCAGCAGCTGCTCCGTTGGTTGCA
GTGATGGGGATGAAAAATACGAGAAGACAAAGTTAGCAGTGGAGATCCAACGTTTTACAAGTTTCATGAA
GAGGATCTCTGTTTGCCAGGGGCAGATTCTGAGGTATTGAGTGGAGTGGAGAGCCTCTCTTCTGCTCTGC
CCCACATTTGAAGTCTCGGAGGTCACAGTCTGCAAGTGGCTGCGGAGGCCAGAGGGCGTTTGAGTTTCAGC
TTATGCCAGCACTGGTCAAGTCAACAGTACTAATCTAGGTCTCACTGTGGAATTTGGGACGATTCT
GGTTTACACCTGTAAGCAGAGTTGCTGGCCTCAAATCAGCAGATGCCCATGGAAGAATTCTGTGTATA
CAAGAAGACCCAGATGCATTTTTATTTAAGTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001109544



[View online »](#)

Insert Size:	1083 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109544.1, NP_001103014.1</u>
RefSeq Size:	1253 bp
RefSeq ORF:	1083 bp
Locus ID:	689637
Cytogenetics:	1q21