

Product datasheet for **RN210529**

Ucp3 (NM_013167) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Ucp3

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN210529 representing NM_013167
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTTGGACTTCAGCCATCAGAAGTGCCTCCACAACAGTTGTAAGTTCCTGGGGCCGGCACTGCAG
CCTGTTTTGCTGATCTCCTCACCTTCCCCCTGGACCCGCCAAGGTCCGCTGCAGATCCAAGGGGAGAA
CCCAGGAGTGCAGAGCGTGCAGTACCGCGCGTGCTCGGTACCATCCTGACTATGGTGCGCAGAGGGT
CCCCGAGCCCTACAGCGGGCTGGTCTGCTGCGCTGCATCGCCAGATGAGTTTGCCTCCATTGCAATTG
GCCTCTACGACTCTGTCAAGCAGTTCTACACCCCAAGGAACGGACCACTCCAGCGTCGCCATCAGGAT
TCTGGCAGGCTGCACACGGGAGCCATGGCAGTGACCTGTGCTCAACCCACGGATGTGGTGAAGGTCCGA
TTTCAAGCCATGATACGCCTGGAACTGGAGGCGAGAGGAAATACAGAGGGAATATGGATGCCTACAGAA
CCATCGCCAGGGAAGAAGGAGTCAGGGGACTGTGAAAGGGACTTGGCCCAACATCACAAGAAACGCCAT
TGTCATTGTGCTGAGATGGTGACCTATGACATCATCAAGGAGAAGCTGCTGGACTCTCACCTGTTCACT
GACAACTTCCCTGTCACTTTGTCTCTGCCTTTGGAGCTGGTTTCTGTGCCACGGTGGTGGCCTCCCCAG
TGGATGTGGTAAAGACCCGATACATGAACGCTCCCCAGGCAGGTACCGAAGCCCCCTACACTGTATGCT
GAGGATGGTGGCCAGGAGGGCCCCACAGCCTTCTACAAAGGATTATGCCCTCCTTTCTGCGTCTGGGA
TCCTGGAACGTGATGATGTTTGTAACTACGAGCAGCTGAAACGGGCTTGATGAAAGTCCAGGTACTGC
GGGAATCTCCATTTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_013167

Insert Size: 927 bp



OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_013167.2 , NP_037299.1
RefSeq Size:	2395 bp
RefSeq ORF:	927 bp
Locus ID:	25708
UniProt ID:	P56499
Cytogenetics:	1q32
Gene Summary:	catalyzes a proton leak across the inner mitochondrial membrane to uncouple oxidative phosphorylation; may reduce the concentration of reactive oxygen species in mitochondria [RGD, Feb 2006]