

Product datasheet for **RN210633**

Kcns2 (NM_023966) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kcns2 (NM_023966) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Kcns2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

RefSeq: [NM_023966.1](#), [NP_076456.1](#)
RefSeq Size: 1660 bp
RefSeq ORF: 1434 bp
Locus ID: 66022
UniProt ID: [Q9ER26](#)
Cytogenetics: 7q22
Gene Summary: human homolog is a alpha subunit of a delayed-rectifier potassium channel [RGD, Feb 2006]