

Product datasheet for **RN210643**

Wnt1 (NM_001105714) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Wnt1 (NM_001105714) Rat Untagged Clone
Tag: Tag Free
Symbol: Wnt1
Synonyms: Int1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN210643 representing NM_001105714
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGGCTCTGGGCGCTGCTGCCAGCTGGGTTTCTGCTACGTTGCTACTGGCACTGACCGCTCTGCCCG
CAGCCCTGGGCGCAACAGTAGTGGCCGATGGTGGGCATCGTGAACATAGCCTCCTCCAGAACCTGCT
GACAGATTCCAAGAGTCTACAGCTGGTGTAGAGCCAGCCTGCAGCTACTGAGCCGCAAGCAGCGACGA
CTGATCCGACAGAACCCGGGGATCCTGCACAGCGTGAGCGGAGGGCTGCAGAGCGCGGTTTCGAGAGTGCA
AGTGGCAATCCGAAACCGCCGCTGGAACCTGCCCACTGCTCCCGGGCCCCACCTCTTCGCAAGATCGT
CAACCGAGGCTGCCGAGAAACAGCGTTTATCTTCGCAATCACCTCCGCGGGGTACACATTCTGTGGCG
CGCTCCTGCTCCGAAGGCTCCATCGAATCCTGCACCTGCGACTACAGGCGGCGGGTCTGGGGGCCCCG
ACTGGCACTGGGGGGTTCAGCGACAACATCGATTCGGTCGCCTCTTCGGCCGAGAGTTTCGTGGACTC
TGGGGAGAAGGGGCGGGACCTACGCTTCCTCATGAACCTCACAATAACGAGGCAGGGCGCACGACCGTG
TTCTCTGAGATGCGCCAAGAGTGCAAATGCCACGGGATGTCGGTTCCTGCACCGTGCGCACGTGTTGGA
TGCGGCTGCCCACTGCGCGCTGTGGGCGAGCTGCTGCGGACCGCTTCGACGGCGCTCCCGCGTCTCT
TTACGGCAACCGAGGCAACAACCGCGCTTCGCGGGCGGAGCTGCTGCGCCTGGAGCCCCAAGACCCCTGCG
CACAAGCCGCCCTCCCCTCACGACCTGGTCTACTTCGAGAAATCGCCCAACTTCTGCACGTACAGTGGAC
GCCTGGGCACAGCTGGCACAGCTGGGCGAGCTTCAACAGCTCATCTCCAGCGCTGGACGGCTGCGAGCT
GCTGTGCTGTGGCCGAGGCCATCGCACGCGCACGAGCGCTCACCGAGCGCTGCAACTGCACCTCCAC
TGGTGCTGCCACGTCAGCTGCCGCAACTGCACGCACACGCGCTACTGCACGAGTGTCT**ATGA**

ACGGCTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001105714



Insert Size:	1113 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001105714.1, NP_001099184.1</u>
RefSeq Size:	2360 bp
RefSeq ORF:	1113 bp
Locus ID:	24881
Cytogenetics:	7q36
Gene Summary:	mouse homolog promotes development of midbrain dopaminergic (DA) neurons [RGD, Feb 2006]