

Product datasheet for **RN210688**

Nudt17 (NM_001109340) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Nudt17 (NM_001109340) Rat Untagged Clone

Tag: Tag Free

Symbol: Nudt17

Synonyms: RGD1565469

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN210688 representing NM_001109340
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGCTGCCGCACGGTTGCTGCTGCGCCTTGCCGGGCGCCTGGAGTCTGTGAGCTTCACACAGAGTGTGT
 GTAGCCTCCTGGGCGCCAGGCAGGAACCCGGGCCGTGGCACACGCACTGCAGCTTAGAGCGAGGACAGAT
 GGTCTCTCTAGCACTTCATTCCCGGGAGCTTCGGAGAGACTTCCAATCCAGAGACCCCTTTCTGCCCT
 TTTGGGGCTCTGGACCAGCAGCCTGAGATCTCCACAAATAGAGGTGTGGAGCTGGGTGTGGCAGTCATTC
 TTCAATCCAGTGACCAGACTGTATTGTTGACCCGAAGGGCATGTACTCTCCGTGTTTCCCCAATCTGTG
 GGTCCCTCCAGGCGGGCACATGGAACCTGATGAAGAGATACTGGAGTGTGGGCTTCGAGAGCTGTGGGAA
 GAGTGTGGACTACAGCTGCCCCAGAATCAGTTTTCTTCTGTCTTCTGGGGTTTTGGGAGCTGCCTACC
 CTCCTAGGCTGAGCTGGGGTTTCCCAAATACCATCACCTCATTCTCTATGTCTGGTGATCTCCCGGGA
 GTCACAGGAGCAACTCCAGGCCAGGATCCAAGTAAATCCAAATGAGGTGAATGCCTTTATGTGGCTGGGA
 CCAGATGTAGCAGGGGCTGTGGCTGCCACAGAGGATGGGACGAAGGCACCTGGACTGCTCTCCAGGACC
 TACCACTCTGTCTGCACAACCTGAAGTGAAGGACGATGGAGGAACCCAGCCTCTGGTCTGCCCGTGTC
 CACGCTCATGCGAGTAACCCAGCCACAGCAGAAGAGGACAAAGAAAGGGTCAGCACTGGAACCAAGTTT
 GCCCTCAAGCTCTGGCTGCAACATCTGGGCAGATGA

ACGCGTACGCGGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001109340

Insert Size: 876 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109340.1, NP_001102810.1</u>
RefSeq Size:	1118 bp
RefSeq ORF:	876 bp
Locus ID:	502584
Cytogenetics:	2q34