

## Product datasheet for **RN210763**

### **Pdha2 (NM\_053994) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Pdha2 (NM\_053994) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Pdha2  
**Synonyms:** Pdhal  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN210763 representing NM\_053994  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGC**C

ATGAGGAAATGCTGGCCACTGTGCTGTCTCAAGTATTTTCAGGAATGGTCCAGAAGCCAGCTCTCAGGG  
GACTGCTGTCATCTCTTAAGTTCTCAACGACGCCACCTGTGACATCAAGAAATGCGACTGTACCTGTT  
GGAACAGGGCCCCCACCTCCACTGTGCTCACCAGAGAGGAGCCCTGAAATACTACCGAACATGCAG  
GTAAATTCGGCGCATGGAGTTGAAGGCCGACCAGCTGTATAAGCAGAAATTCATCCGTGGTTTCTGTACC  
TGTGTGATGGGCAGGAAGCCTGCAATGTGGGGCTGGAGGCAGGGATAAATCCCACAGACCACATCATCAC  
GTCCTACAGAGCTCATGGCTTGTGCTACACGAGGGGACTGTCAGTGAAGTCCATTCTCGCAGAGCTGACT  
GGACGCAAAGGAGGCTGTGCTAAAGGCAAGGGAGGCTCTATGCACATGTACGCCAAGAATTCTATGGGG  
GCAATGGCATTGTTGGGGCACAGGTACCCTGGGAGCTGGTGTGGCTCTCGCCTGTAATACTGAAGAA  
TGGTCAGATCTGCTTGGCTTTGTATGGCGATGGTGGCCTAACCAAGGGCAGGATTTGAAGCTTACAAT  
ATGTCAGCCTTGTGAAATACCTGTGTTTCATCTGTGAGAAATACCGCTATGGAATGGGAACAGCTA  
TTGAGAGATCAGCAGCCAGTACTGATTACCACAAGAAAGGCTTTGTTATCCCTGGACTGAGGGTGAATGG  
GATGGATATTCTAGTGTTCGTGAGGCGACCAAGTTTGCAGCTGATCACTGCAGATCCGGAAGGGGCC  
ATTGTGATGGAGTTGCAGACATACCGTTATCATGGACACAGTATGAGCGATCCAGGGATCAGTTATCGTA  
CACGAGAAGAAGTTCAGAACGTGAGGAGTAAGAGTGATCCAATAATGCTTCTCCGAGAGAGAATGATAAG  
CAACAACCTCAGCAGTGTGAAGAATTGAAAGAAATTGATGCAGATGTGAAGAAAGAAGTCGAGGAAGCA  
GCTCAGTTTGCTACAACCTGATCCAGAACCCTCTGGAAGATTTAGCCAATTACCTCTACCACAAAATC  
CACCATTTGAAGTCCGTGGTGCACATAAATGGCTGAAGTTAAGTCTGTCAAGTAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI



[View online »](#)

<b>ACCN:</b>	NM_053994
<b>Insert Size:</b>	1176 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_053994.2</a></u> , <u><a href="#">NP_446446.1</a></u>
<b>RefSeq Size:</b>	1574 bp
<b>RefSeq ORF:</b>	1176 bp
<b>Locus ID:</b>	117098
<b>UniProt ID:</b>	<u><a href="#">Q06437</a></u>
<b>Cytogenetics:</b>	2q44
<b>Gene Summary:</b>	mouse homolog codes for one of two E1 alpha subunits for pyruvate dehydrogenase (Pdh) [RGD, Feb 2006]