

## Product datasheet for **RN210899**

### **Tbc1d21 (NM\_001014210) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tbc1d21 (NM\_001014210) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Tbc1d21  
**Synonyms:** RGD1310243  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN210899 representing NM\_001014210  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGACCACCTTTCTCCAGAAAACAGCCTCTCTGCCAGGCGGTACGCCACCTTCATTCTGGAGAAGAAAA  
ACCCCCCATTGACAAGGCAGAATGGGACAGCTTCTTCGATGAGAACGGTCACTTGGCCAAGTCACGAGA  
CTTCATTTGTGTTAACATCCTGGAAAGGGGACTACATCCCCTCGTGAGGGCAGACGCCTGGAAGTTCCTT  
ACAGGGTACTACTCATGGCAGAGTTCGCCGGATGAGCGGCTCATGGTGGACAGCAACAGAAGGAGGA  
ATAACTCCTTGTGCCAGATGTATGAGAAAAATCCAGCCCCTGCTGGAAAAATTTACATGGGAACTTCACGGA  
GACTCGGAACAACATCGTGTACGATATTCAGAAGCTGTATGATAAAGATCCCCTGGGCAATGTGCTTGTA  
GACAAGAAGAACTGGAGAAGACCCTGCTGCTGAGTTATGTCTGCAATACAAAGGCTGAGTACCAGCGAG  
GCTTCCACGAGATGGTATGCTTTTCCAGCTGATGGTAGAGCAGCACCATGAGACCTTCTGGCTATTCCA  
GTTTTTCTGCAGAAAACAGAGCACAGTTGTGTCATCAACATCGGTGTAGGCAAGAACCTGGACATGCTC  
AACAACTGATCACTTTACTGGACCCTCAGTTTGTGCTGACCACCTCAAAGGAAAGGGCTCGGGAGCTGTGC  
AGTCCCCTTTCCATGGTCTGCCTTTGCTTCCAGCGTGCCTTCAAATCCTTCGATGATGCTGGAGGCT  
CTGGGAGGCTCTGCTGACAGGGAAGCCCTGCAGGAACCTCCAGGTGCTGGTGGCCTACAGCATGCTACAG  
ATGGTGCAGGAGCAGGCGCTGCTGGAGTCTATGAGCGGTGACGCCATTCTCATGGCCTGCAACAAGCTTA  
TCGACCTTGATGCGGATGAGCTGATCTCTGCTGCCTGCGTGGTTTATTCTGAGCTCATACAAAAGGAGGT  
TCCTCAGGCATTAAGGATTTCTTTCT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001014210



[View online »](#)

<b>Insert Size:</b>	1011 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001014210.1</a></u> , <u><a href="#">NP_001014232.1</a></u>
<b>RefSeq Size:</b>	1249 bp
<b>RefSeq ORF:</b>	1011 bp
<b>Locus ID:</b>	363072
<b>Cytogenetics:</b>	8q24