

Product datasheet for **RN210909**

Fam71f1 (NM_001024324) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam71f1 (NM_001024324) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Fam71f1
Synonyms:	Tdrns18l1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN210909 representing NM_001024324 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATGACATCAGTTCCACCTAGAAAGACTTGGTGGAAATCGAAGAAGACAGTAAAAGTGAGAAGATCTT
ATTCAACCTTCCCTTCTTGAATGCCTGGGAAAAATTCAGGGGCCTCTGCCTGTGGATGGGAAACGAA
CCCTGGAGTGGGCTGGGTGTGGAAGAAGGACTACTCTGCCAGATGGTTCATTCTCCAGAATTCAATCTG
TTTCCTAACTCTGTGGTGTGAAAGCAACTTTGTCCAGGTGAGAAGGAGCAGAGACTGGAAAGAAATCT
ATAAGGCTTCCAACACCATGGCCCTTGGGGTGACCTCCTCCGTACCTTGCTTCCCCTTCCGAACATCCT
CCTTATGGCCAGAGTCATATGGCATCAGGGGCAGAGCCAGACATGGAACAGACCGTCCACAGTGCCAAGC
ATCAACCTGAAGAGCATTCTCCCTTTGAAGTTTGTGGAGCTCCAGATCTGGGACCACCATGAACGCATCC
TACGGTTGAGGACCGTCACCGAGAAAACTACTTCCTAAAGCTCCACCCTGACCATCCTGAGACTGTCTT
CCGCTTTTGGATCCGCCTGGTTCAAATTCTGCATAAGGGTCTGTCCATCACCACCAAGGACCCTACCATC
CTCGTCACTCACTGCCTGGTACCAAGAGCCTTTCAGCCAGGTGGGAAGTCTGAGTTAGTGCAGAAGA
ATTCTAAAGGCCTCCAGCCAGCGAGAGCCTCACACATCTGATGGCCAGGGGAGAGCGAGGCCCTCTC
CCAGATCTTTCCGACTTGCATCAGCAGAAACAGTACAGTTCAGCAGGAGTGAAAAGGTGCAATCAAC
AAGACCAGCTCAGAGAAAGCTACTCCCTGTGAAGACAGTATCCCTGTACCTGCGATCTCAATTGGAGGG
ATGCGTTTCATGTTTCGGAGAGTGGGAAAGAGAGAACCCCTCTGGGCCACAGCCCACTCGCTCCTCAGTAC
CCTGGCTGCCTCCAGCAGGCCACGGCTGTCCCTAATTGGAGGAACTCTAT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001024324
Insert Size:	1035 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001024324.1</u> , <u>NP_001019495.1</u>
RefSeq Size:	1600 bp
RefSeq ORF:	1035 bp
Locus ID:	500061
UniProt ID:	<u>Q68FV5</u>
Cytogenetics:	4q22