

## Product datasheet for **RN210965**

### **PPP1ca (NM\_031527) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ppp1ca (NM\_031527) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Ppp1ca  
**Synonyms:** PP1alpha  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN210965 representing NM\_031527  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGTCCGACAGCGAGAAGCTCAACCTGGATTCCATCATCGGGCGCCTGCTGGAAGTGCAGGGCTCACGGC  
CTGGAAGAATGTGCAGCTGACAGAGAACGAGATCCGTGGTCTTTGCCTCAAATCCCGGGAGATTTCT  
GAGCCAGCCTATTCTTCTGGAGCTTGAGGCGCCTCTCAAGATCTGTGGTACATCCATGGCCAGTACTAT  
GACCTTCTACGGCTGTTTCGAGTATGGTGGCTTCCCTCCAGAGAGCAACTACCTCTTCTTGGGGACTATG  
TGGATCGGGAAAAGCAGTCTTTGGAGACCATCTGCTTGTGTCTGGCCTATAAGATCAAATACCCGGAGAA  
TTTCTTCTACTACGTGGGAACCATGAGTGTGCCAGCATCAACCGCATCTATGGCTTCTATGATGAATGC  
AAGAGACGATACAACATCAAAGTGTGGAAGACTTTCACCGACTGCTTCAACTGCCTGCCATTGCAGCCA  
TTGTAGATGAGAAGATCTTCTGCTGCCACGGAGGCTGTCTCCAGACTTGCAATCCATGGAGCAGATTAG  
ACGTATTATGCGGCCACGGACGTGCCTGACCAGGGCCTGTTGTGTGATCTGCTGTGGTCTGACCCTGAC  
AAGGATGTTCAAGGCTGGGGCAGAAATGACCGTGGAGTCTCCTTACCTTTGGGGCCGAGTTGTAGCCA  
AGTTCCTGCACAAGCATGATTTGGACCTCATCTGCAGAGCACATCAGGTTGTAGAAGATGGCTATGAGTT  
CTTTGCCAAGAGGCAGCTGGTGACACTTCTCAGCTCCCAACTACTGTGGCGAGTTTGACAACGCTGGC  
GCCATGATGAGTGTGGACGAGACTCATGTGTTCTTCCAGATCCTCAAGCCCGCTGATAAGAATAAGG  
GGAAGTATGGGCAGTTCAGTGGCCTGAACCCGGAGGCGCTCCCATCACTCCACCCCGCAATTCTGCCAA  
AGCCAAGAAATAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_031527



[View online »](#)

<b>Insert Size:</b>	993 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_031527.1</a> , <a href="#">NP_113715.1</a>
<b>RefSeq Size:</b>	1392 bp
<b>RefSeq ORF:</b>	993 bp
<b>Locus ID:</b>	24668
<b>UniProt ID:</b>	<a href="#">P62138</a>
<b>Cytogenetics:</b>	1q43
<b>Gene Summary:</b>	plays a role in regulation of synaptic transmission [RGD, Feb 2006]