

## Product datasheet for **RN211030**

### **Sgca (NM\_001107039) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Sgca (NM\_001107039) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Sgca  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN211030 representing NM\_001107039  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGTAGCAGCCGCAACTTGGATTCTCTCCTTGCAAGTCTCCTGGCAGGCTGAGGGACAGCAAGGCC  
AGCAGACAACTTTACACCACTTGTGGTCTGTGTTTGTGCATCCTTTGGAACATGCTACCTTCTGCG  
CCTTCCAGAACACATTGCGGTGCCACCACTGTCCGACTCACCTACCATGCTCACCTCCAGGACATCCA  
GACCTGCCAGGTGGCTGCACTACACACAACGCAGTCCCTATAACCCTGGCTTCCCTATGGCACTCCCA  
CTCCGGAAGATCGTGGGCACCAAGTCATTGAGGTCACAGCCTACAATCGTGACAGCTTTGACACCACTAG  
ACAGAGGCTGCTGCTGCTGATTGGGGACCCGAAGGTCCCAGGTTACCATACCAAGCTGAATTCCTGGTG  
CGTAGCCATGATGTGGAGGAGTCTGCCCTCCACACCTGCCAACCCTTCTCACCGCCTTGGGGGAC  
TCTGGGAGCCAGGGGAGCTCCAGTTGCTCAACATCACCTCTGCCTGGACCGAGGAGGCCGAGTCCCTCT  
TCCTATTGAGGGCCGGAAGGAAGGGTTTACATTAAGGTGGGCTCTGCCACACCTTTTCCACCTGCCTG  
AAGATGGTGGCATCACCTGACAGCTATGCCGTTGTGCCAGGGACGCCTCCACTATTGCTCTGCTATG  
ACACCTTGGCACCTCACTCCGTGTCGACTGGTGAATGTATCTCTGGTGGACAAGTCAGTGCCAGAGCC  
CCTGGATGAGGTGCCTACTCCTGGCGATGGATCTTGGAGCACGACCCATTCTTCTGCCCGCCACTGAG  
GCCACGGGCCGAGACTTCTTAACAGATGCCTTGGTGACCTCTTAGTGCCCGTGCTGGTGGCTCTGCTGC  
TCAACCTGTTGCTGGTTTACATTATGTGCTGTCCGGCTGAAGGACGGCTGAAGAGAGACATGGCCACCTC  
TGACATCCAGATGTTTACCACCTGCTCCATCCACGGGAATACGGAAGAGCTTCGGCAGATGGCAGCCAGC  
CGCGAGGTGCCCCGGCTCTATCCACCTTGCCCATGTTCAACGTGCGTACAGGAGAGCGGTTACCTCCCC  
GAGTAGACAGCGCACAGATGCCCTCATCTGGACCAGCA**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001107039



<b>Insert Size:</b>	1164 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001107039.1</a></u> , <u><a href="#">NP_001100509.1</a></u>
<b>RefSeq Size:</b>	1693 bp
<b>RefSeq ORF:</b>	1164 bp
<b>Locus ID:</b>	303468
<b>Cytogenetics:</b>	10q26