

Product datasheet for RN211048

Gna11 (NM_031033) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gna11 (NM_031033) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Gna11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN211048 representing NM_031033 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACTCTGGAGTCCATGATAGCGTGTTCCTGAGCGACGAGGTGAAGGAGTCAAGCGCATCAACGCGG
AGATCGAGAAACAGCTGCGGCGAGACAAGCGGACGCCCGGCGGAGCTCAAGCTGCTGCTTCTCGGCAC
TGGCGAGAGCGGGAAGAGCACCTTCATCAAGCAGATGCGCATCATCCATGGCGCCGGCTACTCGGAGGAG
GACAAACGCGGCTTCACAAAGCTGGTGTACCAGAACATCTTACCGCCATGCAGGCTGTGGTGCGCGCTA
TGGACACGCTCAAGATCCGCTACAAGTATGAGCAGAACAAGCCAATGCACTCCTGATCCGGGAGTTGA
TGTGGAGAAGGTCACGACTTTTGAGCATCAGTATGTGAACGCCATCAAGACACTGTGGAGTGACCCCTGGT
GTCCAGGAGTGCTACGATCGCAGGCGAGAGTTCAGCTCTCTGACTCGGCTAAGTACTACTTGACGGACG
TGGACCGCATCGCTACAGTAGGCTACCTGCCACCCAGCAGGATGTGCTGCGGGTACGTGTGCCACAAC
TGGCATCATCGAGTACCCGTTTGACCTGGAGAACATCATCTTCAGGATGGTGGATGTGGCGGGCAGAGG
TCAGAGCGCAGGAAGTGGATCCATTGCTTTGAGAATGTCACGTCCATCATGTTCTTGGTGGCATTAAAGC
AGTATGACCAGGTCCTGGTGGAGTCGGACAACGAGAACCAGGAGGAGCAAGGCCCTGTTCCGCAC
AATCATCACCTACCCCTGGTTCAGCACTCGTCTGTGATCCTCTCTCAACAAGAAGGACCTTCTGGAA
GACAAGATCCTGCACTCACACTTGGTCGACTACTTCCCGAGTTCGATGGCCACAGAGGGACGCACAGG
CCGCCCGGAGTTCATCCTGAAGATGTTTGTGGATCTGAACCCGACAGCGACAAGATCATCTACTCGCA
CTTCACCTGCGCCACGGACCCGAAAACATCCGCTTCGTGTTGTGCCGTGAAGGACACCATCCTCGAC
CTGAACCTGAAGGAGTACAACCTGGTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_031033



Insert Size:	1080 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_031033.1 , NP_112295.1
RefSeq Size:	1080 bp
RefSeq ORF:	1080 bp
Locus ID:	81662
UniProt ID:	Q9JID2
Cytogenetics:	7q11
Gene Summary:	may be involved in mediating second messenger signals from the endothelin receptor in adipocytes [RGD, Feb 2006]