

Product datasheet for **RN211080**

Zbtb49 (NM_001107223) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zbtb49 (NM_001107223) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Zbtb49
Synonyms:	Zfp509; Znf509
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >RN211080 representing NM_001107223
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGC**

ATGGACCCTGTTGCTATTCATAGCTGCCACCTCCTCCAGCAGCTGCGAGAGCAGCGGATCCAGGGCCTGC
 TTTGTGACTGCATGCTGGTGGTGAGGGCGTGTGCTTCAAAGCACACAAGAATGTGCTGGCAGCCTTCAG
 CCAGTATTTTCAGGAGTCTTTTTCAGAACTCTTCCAGCCAGAAGAATGACATCTTTACCTGGACGTTAAA
 AACGTCAGTGGCATCGGGCAGATCCTGGACTTCATGTACACGTCGCGCCTCGACCTTAACAGGACAACA
 TTCAGGTGATGCTGGACACAGCGCAGTGTGCAAGTGCAGAACGTTCTGAATCTGTGCCACACGTTTCT
 GAAGTCGGCCTCCGCGGCACCGCTGCCAGGCTTGCCTTGTGCTGGCGTTTCTCCCTGCAGGGTGTGAGT
 CTGGACGGCACCTGCGCGGCGGAGTGAGCACTACCCCTCACTCACTGCAGGACTGCCCTATGGAAG
 GCCCGCGGGCTAAAGACGTGTCCAGAGGGAATCCCGTGTCTCCACAGCAGATTTCCGGTCGTCCACAGG
 AGAAGGGTCAAAACAGCTGCCCAAGTGGCAGCTGCCAGAGGGCGCCCGCAAACAGCCCAACTGCTAT
 TACAAACTCAGAACACTGTACAGCAAGCAGTACTACAAGCAAAGTGCCTGCCCCAGTCCGGTGCCTGCCA
 CAGAGCAGCCGCCCACTCCCGTGCCTCCACAGACCTCGTGCAGCAGACTCACAGCCGCTGTGGAGGG
 CCGTGCAGGGGTCCCGGAGTCTCCAGAGCACTTGGCTTCACTGTTGTGGCTCCACCGGTGAGGAGCTCT
 GGCATCGACTCAGAGGCAGACACCTTTCTCAGCCACCTGCCAAGCAGATGAGGCTGAAGAAGGCCATGC
 ACCTGAAGAACTAACTTTCTTAAGTCACAGCAGTCACTGAGTGCACCTCACATCCTCAGTCAGATAA
 TGTGTTAGTCAGGGGGAGGAATCTGCTGCTAAGGAAGATGCAGGGGAAAGAGCCGGAAGCCAGACAAC
 GAAGACAAAGGGAGGGAAGAAGTGGTCCAGAGAGCAGCCGCGAGGTGGAGATGCCAGGAGCCCCAGCTT
 CCTGGGAAGACCGTCCAGGCCCTCCAGTCCAGAAACAGTATGCGTGTGAATTGTGTGGGAAGCCTTT
 TAAACACCCAAGCAATCTGGAGCTGCACAAACGCTCTCATACAGGCAGGCAACCTTCAGACTCACTTACG
 CCGGCACTCTGGGAGAAGCCATACATCTG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001107223

Insert Size: 1293 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001107223.1, NP_001100693.1

RefSeq Size:	2685 bp
RefSeq ORF:	1293 bp
Locus ID:	305428
Cytogenetics:	14q21