

Product datasheet for RN211122

C1d (NM_001106021) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C1d (NM_001106021) Rat Untagged Clone
Tag:	Tag Free
Symbol:	C1d
Synonyms:	RGD1560600
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN211122 representing NM_001106021 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCAGGTGAAGAAATGAATGAAGACTATCCCGTAGAAATTCACGATTCTTTAACAGCTCTGGAGAGCT
CCCTGGGGCGGTGGATGACATGCTGAAGACCATGATGTCTGTTTCTAGAAATGAGTTGTTGCAGAAAGTT
GGACCCACTGGAACAAGCAAAGGTGGATTTAGTTTCTGCATATACATTAAATTCAATGTTTGGGTTTAC
TTGGCAACTCAAGGAGTTAATCCCAAGGAGCATCCAGTGAAGCAGGAGTTGGAAGAATCCGAGTCTACA
TGAACAGAGTTAAAGAAATAACAGACAAGAAGAAGGCTGCCAAGCTGGACAGAGGGCCGCGTCGAGATT
TGTCAAAAATGCACTCTGGGAACCCAAACCAAAAAACACCAAAAAGTGGCTAATAAAGGAAAAAGTAA
CACT**AA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001106021
Insert Size:	426 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106021.1, NP_001099491.1</u>
RefSeq Size:	3129 bp
RefSeq ORF:	426 bp
Locus ID:	289810
Cytogenetics:	14q22