

Product datasheet for **RN211190**

Lrwd1 (NM_001013981) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lrwd1 (NM_001013981) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Lrwd1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN211190 representing NM_001013981 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCTCCTCTCACGCCGAGCTGCTCCTGCAACGTGGGCGACCCAAGACCGACAGGCTGGGGAAGATTC
AGAGTCTGAATCTGTGAGGGCTGGAGCTGCTCTCAGAACACTTGGATCCCAATCTCCTGGGACGCCTGAA
GAAGCTGAAGGAACCTGACCTCTCTAACAACTTGTGGAACATTGCCTGCTAACCTAGGCCTGTCCCAT
CTGCGAATCCTTCGCTGTACCAACAACAGCTAGGGGACGTCACCACCTTGACACAGTTCCTTCCTGCCC
AGGAGCTCAGTCTGGAAGGCAACCTTTCTGACGGTCAGTGACAACCTCAAAGTCTCCTTCTCCTGCC
CAAGCTGCGTAAGGTCAATGGCAAGGACACGGCCTCCACCTGCTCACAGGTCACAGCCCTCTGGCAGAAG
TTCATAGCCACTGTGAACCTGAGGAGGAAGCTAACAAAGTGCAGGCCGACTTCATGAGGTCCGCAGTCA
GGGATGTGCGTTACGGGCCGAGTCTCTCTCGAGTTCACGCAGTGGCGGGTGGGATGATCGCTGAGGA
ACTGGTGGCTCAGGTGGGGCCAGGCCACGAAGCTAATGCATCAGTGAACATCCCCAGGCAGCAGTT
GGCTCCAAGCTCAGGGCCAGAAAGGTAACAAAACGTCCAGATGATGACGAAATCAGCCTTCTCCCACCA
AGCGAGTGGGGCCTCCCCACGAGCCAGACTGAGGACAGTCTATGGATGCTGCTGGAGGCCAGGCCGC
CCTACACCTGGAGCCCTGCACCTCTGTCAGTGCCACAGCAGAAACAACAGCCCTAAGGACCTGGAGACC
CAGCTGTGGGCTGTGCCTTTGAGCCGGCCGGGAGGAGGACACTCGGGGGCCACATCCCAGACCGTGG
CCACATGTGGCGGAGAAGCTGTTTGCCTGATAGACTGCCAGACAGGCCTCGTGCTCCATAAATAAAGGT
TCCGGGCGAGGAGTTCTTCTCAGTGGCCTGGACAGCCCTGACGGTGGCCACCCAGGCAGGCCACAAGAAG
CGCTGGAACATGCTGGCAGCTGCAGGCCTTCGGGGCATGGTGGGCTGCTCCATGTGCGTGCGGGCTTCT
GCTGCAGCGTCATCAGGGCCCATAGAAGGCCATCGCCACCCTCTGCTTCAGCCCCACCATGAGACTCA
CCTCTTCAGTAAGACCCTCAGCTCTGGGCTCTCCAGCACTATAGTCTTCTGCCTCG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



View online »

ACCN:	NM_001013981
Insert Size:	1251 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001013981.1, NP_001014003.1</u>
RefSeq Size:	2564 bp
RefSeq ORF:	1251 bp
Locus ID:	304396
Cytogenetics:	12q12