

Product datasheet for RN211220

Htr6 (NM_024365) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Htr6 (NM_024365) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Htr6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN211220 representing NM_024365 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTTCCAGAGCCAGGCCCTGTCAACAGTAGCACCCAGCCTGGGGTCCCGGGCCACCGCCTGCTCCGG
GGGGCAGCGGCTGGGTGGCTGCCGCGTGTGCGTGGTCATCGTGTGACAGCAGCCGCAATTTCGCTGCT
GATCGTGCTCATTGACGCAGCCCGCGCTGCCAACACGTCTAACTTCTTTCTGGTGTGCTCTTCAGC
TCGGACTTGATGGTGGGGTTGGTGGTATGCCCCAGCCATGCTGAACGCGCTGTATGGGCGCTGGGTGT
TAGCTCGAGGCCCTGTCTGCTTTGGACTGCCTTCGACGTGATGTGCTGCAGCGCCTCCATCTCAACCT
CTGCCTCATCAGCCTGGACCGCTACCTGCTCATCCTCTCGCCGCTGCGCTACAAGCTGCGCATGACAGCC
CCGCGAGCCCTGGCGCTCATCCTGGGTGCCTGGAGCCTCGCGGCGCTTGCTCCTTCTACCCCTCTTGC
TGGGCTGGCAGCAACTGGGCAAAGCTCGAACACCTGCCCTGGCCAGTGGCCCTATTGGCCAGCCTGCC
TTTTGTCTCGTGGCGTCCGGCGTCACTTTTTCTGCTTTCGGGTGCCATCTGCTTCACTACTGCAGG
ATCCTTCTGGCTGCCCGAAGCAGGCGGTGCAAGTGGCTCGCTCACCACGGGCACGGCTGGCCAGGCT
TGGAAACCTTGCAGGTGCCAGGACACCAGCCAGGGATGGAGTCCGCTGACAGTAGGCGTCTGGCCAC
CAAGCATAGCAGGAAGGCCTTGAAGCCAGCCTGACCTGGGCATCCTGCTGGGAATGTTCTTTGTCAAC
TGCTGCCCTTCTTTGTGGCCAACATAGCTCAGGCCGTGTGTGACTGCATCTCCCAGGCCCTCTTCGATG
TCCTCACATGGCTGGGGTACTGTAAATAGACCATGAACCTATCATCTACCCGCTCTTATGCGGGACTT
CAAGAGGGCCCTGGGCGAGTTCCTGCCATGCGTCCACTGTCCCCGGAGCACCGGGCCAGCCCTGCCTCC
CCCTCCATGTGGACCTCTCACAGCGGTGCCAGACCAGGCCCTCAGCCTGCAGCAGGTGCTGCCTCTGCCTC
TGCCGCCAAACTCAGATTCAGACTCCGCTTCAGGGGGCACCTCGGGCCTGCAGCTCACAGCCAGCTTCT
GCTGCCTGGAGAGGCCACCCGGGACCCCGCCACCCACCAGGGCCACCACTGTGGTCAACTTCTTTGTG
ACAGACTCTGTGGAGCCTGAGATACGGCCGCATCCACTCAGTTCCCCCGTGAAC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_024365
Insert Size:	1317 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_024365.2 , NP_077341.2
RefSeq Size:	1925 bp
RefSeq ORF:	1317 bp
Locus ID:	64354
Cytogenetics:	5q36
Gene Summary:	G-protein coupled receptor for serotonin [RGD, Feb 2006]