

Product datasheet for RN211230

Eif3l (NM_001034134) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Eif3l (NM_001034134) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Eif3l
Synonyms:	Eif3s6ip
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN211230 representing NM_001034134 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGCCCTGCCACATTGTTTCCTTCTGCAACAGTTCAGTCATTCACTCAGTCAGTATCGCTGCAAGACGGCCA
 AGAAGTCAGAGGAAGAGATCGACTTCCTGAGATCCAACCCTAAGATCTGGAACGTGCACAGCGTCTCTCA
 CGTGCTGCATTCTCTGGTGGACAAATCCAACATCAACCGGCAGCTGGAGGTGTACACCACTGGAGGTGAC
 CCCGAAAGTGTGGCTGGGAATATGGACGACACTCCCTCTACAAGATGCTTGGCTACTTCAGCCTGGTGG
 GGCTTCTGCGCCTACATTCTCTACTGGGGATTACTACCAGGCTATCAAGGTGCTGGAGAATATTGAGCT
 AAACAAGAAAAGCATGTATTCCCGTGTGCCTGAGTGCCAGGTCAACACCTACTATTATGTTGGTTTTGCA
 TATTTAATGATGCGTCGATACCAAGGATGCTATCCGCGTCTTTGCCAACATCCTCCTGTACATCCAGAGGA
 CCAAGAGTATGTTCCAGAGGACCACATACAAGTATGAAATGATTAACAAGCAGAATGAGCAGATGCATGC
 TCTGTTGGCCATTGCTCTACCATGTACCCATGAGGATTGATGAGAGCATACACCTCCAGCTGCGGGAG
 AAATATGGCGACAAGATGCTGCGCATGCAGAAGGGCGACCCCAAGTCTATGAGGAACTTTTAGCTATG
 CCTGCCCAAGTTCCTGTGCGCTGTGGTGCCTAACTACGACAACGTGCATCCTAACTACCACAAAGAGCC
 CTTCTGCAGCAGCTGAAGGTGTTTTCTGATGAAGTGCAGCAGCAGGCCAGCTCTCCACCATCCGCAGC
 TTCCTCAAGCTCTACACCACCATGCCTGTGGCCAAGCTGGCTGGCTTCTTGACCTCACAGAGCAGGAGT
 TCCGATCCAAGTGTCTTCTTCAAGCACAAGATGAAGAACCTGGTGTGGACCAAGTGGCATTCTCTGCCCT
 AGATGGCGAATTCCAGTCGGCCTCGGAGGTGGACTTCTACATTGATAAGGACATGATCCACATCGCAGAC
 ACCAAAGTTGCCGACGCTATGGGGATTCTTCATCCGGCAGATCCACAAGTTTGAAGAGCTTAATCGAA
 CCTGAAGAAGATGGGACAGAGGCC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001034134
Insert Size:	1149 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001034134.1, NP_001029306.1</u>
RefSeq Size:	1514 bp
RefSeq ORF:	1149 bp
Locus ID:	300069
Cytogenetics:	7q34