

Product datasheet for RN211235

Mesp1 (NM_001107531) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Mesp1 (NM_001107531) Rat Untagged Clone
 Tag: Tag Free
 Symbol: Mesp1
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Fully Sequenced ORF: >RN211235 representing NM_001107531
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCAGCCCCTGTGCGGCCACTCTCCAAGTCCTGGATCCCAGTCCCCTCGTCCGCAGCCATCGA
 TGCCTCCGATGGGGACAGCTTCTGCTCCCCAGCCTGGTCTCGGACTCGTGGGACGGTCCCAGGCCAG
 CAGCCCTGCACTACCCTGTGCCCGCCCGCCGGCGTGTGGGACCCCGGTAGGCGCGGGACGCATGGC
 AGCCGCTGGGTAGCGGACAGCGGCAGAGCGCCAGCGAACGCGAGAAGCTTCGTATGCGCACGCTCGCCC
 GCGCGCTGCACGAGCTGCGCCGTTCTTGCTCCTTCAGTGGCACCAATCGGCCAGAACCAGCAAGAT
 CGAGACGCTGCGCCTGGCCATCCGCTACATCGGCCACCTGTGCGGTGTGCTGGGCCTCAGCGAGGACAGC
 CTCCGGCAACAGCGACACGCGGTGTACCTCGAGGCTGCCCGCTGTGCCCGACAGCGGCCTGGCGCAGG
 CGCAGTCTCTCGGTCCCATTTAAGCCCGGTTGCTGCAGCGGGGTGCTCTGGGGTTCTTCGCCCGCCTA
 CCCTAGACCCCGGTCGCCGAGAATCGTGGGACCCATCATTCTGTACACCGAAACAGCATCCCTGGAA
 AGGCAGGAGATGGAGCCAGCCCTCGTCTCCGCTCTTCAGCAGCGACATGCTGGCTCTGCTAGAAACCT
 GGACACCTCAGGAGTGGCCGCCAGCTGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
 ACCN: NM_001107531
 Insert Size: 729 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001107531.1, NP_001101001.1</u>
RefSeq Size:	1275 bp
RefSeq ORF:	729 bp
Locus ID:	308766
Cytogenetics:	1q31