

Product datasheet for **RN211248**

Cd209d (NM_001108849) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Cd209d (NM_001108849) Rat Untagged Clone

Tag: Tag Free

Symbol: Cd209d

Synonyms: Clec4m; Dc-signr; RGD1561104

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN211248 representing NM_001108849
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**C

ATGAGTGAGTCCATGGAATCAAACACACAGCAACTGGACAATCCAGAGGATGAAGAGTATTTGATGAGTG
GTTCCAGGTATTCTATTATCAGCTCTCGATTACGAACAAAATTTGGAATCAAAAGTTTGGCAGAATATAC
GAAGCAGAGCCACAACCCCTTGGTCCTACCGCTGCTCTCCTTCTGTTCTTGGCTGGCCTTCTGCTGATC
ATTCTTGTCTAGTCTCAAAGGTACCCAGCTCAGAAGTTCAGGACAAAATCTACCAGGAGCTCATGCAAC
TGAAGACTGAAGTCCATGATGGCTTGTGTCAACCTGTCCGAGGGATTGGACATTCTCAATGGAAGCTG
TTACTTCTTCTCCAAGTCCCAAAGAACTGGCACAACCTCCATCACTGCCTGCAAGGAACTGGGGGCTCAA
TTGGTCATCGTAGAGACGGATGAAGAGCAGACCTTCTGCAGCAGACTTCTAAGACTAGGGGACCAACAT
GGATGGGTCTCTCAGACATGCATAATGAAGCCACGTGGCACTGGGTGGATGGCTCACCTCTGTCACCCGAG
CTTTGCACAGTATTGGAATAGAGGGGAGCCCAACAATGTCGGTGATGAGGATTGTGCAGAGTTCTCTGGG
GACGGCTGGAATGACCTCAGATGCGATACCCGGATTTTCTGGATCTGTAAGAAAGTTTCAACCTCATCAT
GTACCATCAAA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001108849

Insert Size: 714 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108849.1, NP_001102319.1</u>
RefSeq Size:	915 bp
RefSeq ORF:	714 bp
Locus ID:	363856
Cytogenetics:	12p12