

Product datasheet for **RN211315**

Peli1 (NM_001100565) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Peli1 (NM_001100565) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Peli1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN211315 representing NM_001100565 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTTTCTCCTGATCAAGAAAATCATCTTCCAAAGCCCCAGTAAAATATGGCGAACTCATTGTCTTAG
GGTATAATGGGTCTCTCCAAACGGCGACAGAGGAAGGAAAAGTAGGTTTGCTTTGTTAAAAGACC
TAAGGCAAATGGGGTGAAGCCTAGCACCGTGCACATTGCATGCACTCCTCAGGCCGCAAGCAATAAGC
AACAAAGGACCAGCATAGCATATCGTATACTTTATCTCGAGCCCAGACAGTGGTGGTTGAATACTCACG
ACAGCAACACTGATATGTTTCAGATTGGCCGGTCAACTGAAAGTCTATTGATTTGTGGTAACTGACAC
AGTTCTGGAAAGTCAGAGTAATTCGACACGCAGTCCGTACAAAGCACTATATCGAGTTTGCCTGTAGG
ATCATTGTGAACGCAATCCTCCCTTACAGCTCGGATTTATGCTGCAGGATTTGATTTCAAAAAACA
TCTTTCTGGGGAGAAGGCTGCCAAATGGAAGACGTCCGACGGGCAGATGGACGGCCTGACCACCAATGG
AGTTCTGTGATGCATCCACGCGATGGGTTACAGAAGACTCCAAGCCTGGAATATGGAGAGAGATATCC
GTATGTGGGAACGCTTTCAGTCTGCGTGAACCAGATCAGCGCAACAGAGAGGAAAGATGGTGGAAATCG
AAACCAATCAGCTGCAAGATGGCTCCTTAATCGACCTCTGTGGTGAACCCTGCTGTGGCGCACAGCAGA
AGGCCTTTCCACACCCCTACTGTGAAACACTTAGAAGCTCTGAGACAGGAGATCAATGCAGCTCGGCCT
CAGTGCCTGTAGGCTTCAACACTAGCCTTCCCAGCATGAAGAGGAAAGATGTTGTAGATGAGAAGC
AACCATGGGTATATCTGAACTGCGGCCATGTTTCATGGTTACCATAACTGGGAAACAAGGAAGACGTGA
CGGCAAAGATCGTGAATGTCTATGTGTAGGTCTGTTGGTCCCTATGTCCCTCTGTGGCTGGATGTGAA
GCTGGATTTATGTGGACGCCGGCCCTCCACCCATGCCTTTAGCCCCTGTGGGCATGTGTGTTAGAAA
AGACAACGGCTTACTGGTCCCAGATCCCCTTCTCATGGTACGCACACTTTTACGCAGCCTGTCCCTT
CTGTGCACATCAGTTGGCTGGTGAACAAGGCTATATCAGACTTATTTTCAAGGACCTCTAGACTAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001100565
Insert Size:	1257 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001100565.1</u> , <u>NP_001094035.1</u>
RefSeq Size:	2383 bp
RefSeq ORF:	1257 bp
Locus ID:	305549
Cytogenetics:	14q22