

## Product datasheet for RN211377

### Rnf182 (NM\_001109117) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rnf182 (NM\_001109117) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Rnf182  
**Synonyms:** RGD1560399  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN211377 representing NM\_001109117  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCCAGCCAGCCGCCGAAGAGCCTGCAGAGTTTCAGGTCTCAGATGAGCTCGAGTGTAAAGATCTGTT  
ACAACCGGTACAACCTGAAACAGCGGAAGCCCAAGGTCTGGAGTGTGTACAGGGTTTGCCCAATG  
CCTCTACAAGATCATAGACTTTGGGGACTCTCCCAAGGTGTCATCGTCTGCCCGTTCTGCAGGTTTGAG  
ACGTGCCTGCCGGATGACGAAGTTAGCAGTCTGCCCGATGACAATAACATCCTTGTAAATCTGACTTGCG  
GAAGCAAAGGGAAGAAATGCCTGCCCGAGAACCCACCGAGCTGCTGCTTACCCCAAGAGGCTGGCTTC  
CCTTGTACGCCCTTCCCACACATCCTCCAAGTCTTGGTATCACCATCATGGAGGTACAGAGGGAGAGC  
TCCCGTCTCTCAGCTCCACCCCTGTGGTGAATTCTACAGGCCAGCAAGTTTCGACTCTGTACCACAG  
TGTCCACAATTGGACAGTGTGGAAGTGCACCTCCTTGCCTTTTCAGACGTCCATCCGGTACTGGTGTG  
GTTGCTAGGTTTGTGTAAGTCTCAGCTCCTTGCCTTAGGGATCTACTTACTGGTGTCTAAGAAGTCAAC  
CTGGGGGTGGTCTTTGTAGCCTTGTCCCCTTAGCCTTGTATCCTTATGGTATATGGCTTTTGCCAA  
CGTTTTGTCATGAATTTCTGGACTGTATGGCACTTCCTTCT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001109117  
**Insert Size:** 744 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001109117.1</a></u> , <u><a href="#">NP_001102587.1</a></u>
<b>RefSeq Size:</b>	1220 bp
<b>RefSeq ORF:</b>	744 bp
<b>Locus ID:</b>	498726
<b>UniProt ID:</b>	<u><a href="#">D3ZBM4</a></u>
<b>Cytogenetics:</b>	17p12
<b>Gene Summary:</b>	E3 ubiquitin-protein ligase that mediates the ubiquitination of ATP6V0C and targets it to degradation via the ubiquitin-proteasome pathway.[UniProtKB/Swiss-Prot Function]