

Product datasheet for RN211425

Azgp1 (NM_012826) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Azgp1 (NM_012826) Rat Untagged Clone
Tag: Tag Free
Symbol: Azgp1
Synonyms: ZA2GA; Zna2gp
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN211425 representing NM_012826
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGCCTGTCCTGTTGTCCCTGCCTCTCCTTCTGGGTCTGCAGTCCTCAGGAGACTGGGTCTTATT
CTCTGATCTTCTCTACACCGGGTTGTCCAGGCCAGCAAAGGCCTTCTGAGTTTCAAGCCACCGCATT
TCTCAACGATCAGGCCTTCTTCCACTACAACAGCAACAGTGGGAAGGCAGAGCCCGTGGAAACCTGGAGC
CATGTGGAAGGAATGGAGGACTGGGAGAAGGAAAGCCAGTTCAGAGAGCCAGGGAGGAGATCTTCTCTCG
TGACCCTGAAAGACATCATGGACTATTACGAGGACAGTACAGGGTCTCACACCTTTCAGGGAATGTTTGG
TTGTGAGATCACAAATAACAGAAGTAGTGGAGCAGTCTGGAGGTATGCCTACGATGGAGAAGATTTTCATC
GAATTCACAAGGAAATCCCAGCCTGGATCCCCTTAGACCCAGCAGCCGCAAAACCAAGCTGAAGTGGG
AGGCAGAAAAGGTCTATGTGCAGCGAGCCAAGGCATACCTAGAGGAGGAGTGTCTACAATGCTGAAGAA
GTACCTGACCTACAGTAGATCTCACCTGGACCGAACAGATCCTCCACTGTGAAAATCACCAGCCGTGTG
GCCCCAGGAAGAAACAGAATATTCAGATGTCTGGCCTATGACTTCTACCCACAAAGAATAGTCTGCACT
GGAACCAGGCCAGCAAGAAGCTGGCATCTGAACCAGAAAGAGGTGTTTTTCCCAATGGAAATGGCACTTA
CCTCTCCTGGATGGAGGTGGAAGTCCCCCACAGAACAGAGACCCCTTTGTCTGCCACATAGAACACAAG
GGCTATCCCAATCTCTCTCAGTGCAGTGGGATGAGAAAAGTAAAGT**CTGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_012826
Insert Size: 891 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_012826.2 , NP_036958.2
RefSeq Size:	1212 bp
RefSeq ORF:	891 bp
Locus ID:	25294
Cytogenetics:	12q11
Gene Summary:	has similarity to class I major histocompatibility complex (MHC) antigens [RGD, Feb 2006]