

## Product datasheet for RN211427

### Clrn3 (NM\_001014026) Rat Untagged Clone

#### Product data:

|                           |                                                                                    |
|---------------------------|------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                |
| Product Name:             | Clrn3 (NM_001014026) Rat Untagged Clone                                            |
| Tag:                      | Tag Free                                                                           |
| Symbol:                   | Clrn3                                                                              |
| Synonyms:                 | RGD1305043; Tmem12                                                                 |
| Mammalian Cell Selection: | Neomycin                                                                           |
| Vector:                   | pCMV6-Entry (PS100001)                                                             |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                               |
| Fully Sequenced ORF:      | >RN211427 representing NM_001014026<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGCCAACACACAGAAAACGTTGATGTTCTGTCTGGCTTTCTCACAAGCCTGGGCTCAGTTGTTGTCA  
TCTGTTCTATTCTGGGGACACAGGCATGGATCACCAGCAAGATTTCTTCACAGACACTATCTCCAATGG  
TACCATTGCTATCACCTATGGGCTTTTCGTGGTACGAGTAGTCAAGATTTGAACGAAGGCCTTCAAGAG  
CTGGACAAAACTTTGAAGTTTTAGGGATATTGAGCAATCTTCCAGAAGTCCCTGCATGTGGTGGTCA  
TCATACTCCTGATCCTGAGTCTGGCAGCTTCCTTGCTGAGCTCAATGTTACCTTCTACAACAGCATCAG  
CAACCTTACCAGACGTTCTGGGACCTATGGGAGTCTACACCTGGAATGGACTCAGTGCATCCTTTGTG  
TTTCTGACCATGGTTCTGTTTGTGGGAATGTCGACTCCAACCATCTCTCGAAAAGCTGTCCCAGACAC  
TTTATCCAGATGCCATTAACAAGAAAACAACCCACACTTATGGATACTCATTTTGGCTCATCCTGCTTGT  
CATCCTACTGAATATTGTCACCGTGGTCATCATTTTTCTATCAGAAGGCCCGATACCATCAGAAGCAG  
GAGCAGAGGAAGCCAGTGGAGTACGCTCCGAGGGATGGCATTTTATTCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

|                    |              |
|--------------------|--------------|
| Restriction Sites: | SgfI-MluI    |
| ACCN:              | NM_001014026 |
| Insert Size:       | 681 bp       |


[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                    |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                              |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_001014026.1</a></u> , <u><a href="#">NP_001014048.1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq Size:</b>           | 2041 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>RefSeq ORF:</b>            | 681 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Locus ID:</b>              | 309082                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>UniProt ID:</b>            | <u><a href="#">Q6AYR5</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Cytogenetics:</b>          | 1q41                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |