

Product datasheet for **RN211637**

Fam72a (NM_001081451) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Fam72a (NM_001081451) Rat Untagged Clone

Tag: Tag Free

Symbol: Fam72a

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN211637 representing NM_001081451
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCTACCAACAACTGTACTTTCAAAGACCGATGCGTCTCCATCCTGTGTGCAAATTCTGTAACAAG
 TTCTCAGCTCTAGGGGAATGAAGGCTGTCTGTTGGCTGACACTGACATAGACCTTTACTCTACAGACAT
 TCCTCTACCAACACAGTGGATTTCATTGGAAGATGCTATTTACAGGAATCCACAAGTGTAAGTGAAG
 GACATCGCATGTTTAAATGTGGGAACATTGTAGGCTACCATGTGATTGTTCCGTGTAGCTCCTGCCTTC
 TTTCTGCAACAATGGACACTTCTGGATGTTTCACAGCCAGGCGGTTTATGGTATAAACAGACTAGACTC
 GACTGGTGTAACCTTCTACTTTGGGGCAATTTGCCAGAGACAGAAGAATGCACAGATGAAGAAATGCTA
 GAGATTTACGAGAAGAGTATATTCGATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001081451

Insert Size: 450 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001081451.1, NP_001074920.1</u>
RefSeq Size:	450 bp
RefSeq ORF:	450 bp
Locus ID:	681249
UniProt ID:	<u>A1KXW8</u>
Cytogenetics:	13q13
Gene Summary:	May play a role in the regulation of cellular reactive oxygen species metabolism. May participate in cell growth regulation (By similarity).[UniProtKB/Swiss-Prot Function]