

Product datasheet for RN211656

Car11 (NM_175708) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Car11 (NM_175708) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Car11
Synonyms:	CAR-XI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN211656 representing NM_175708 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGGGGTGGGCTCGCCTGAGCGCCCTCAAGCGCTGGTACTCTGGGCCGCACTGGGAGCGGCAGCTC
ACATCGGACCCGCACCTGATCCCGAGGACTGGTGGAGCTACAAGGAGAATCTCCAGGAACTTCGTGCC
AGGGCCTCCCTTCTGGGGCTGGTGAACGCTGCCTGGAGTCTGTGCGCTGTGGGAAGCGGCAGAGCCCT
GTGGATGTGGAGCTGAAGCGGGTCTTTATGACCCTTTTCTGCCCCGCTCAGGCTCAGCACCGGGGGAG
AGAAGCTCGGGGAACCTTATACAACACTGGCCGCATGTCTCCTCCTACCTGCATCCCGCCTGTGGT
CAACGTGTCTGGGGTCCCTTCTATAGCCACCGACTCAGTGAAGTGCAGGCTGCTGTTGGAGCTCGA
GATGGTGTGGCTCCGAACACCAGATCAACCATCAAGGCTTCTCTGCTGAGGTGCAGCTAATCCACTTCA
ACCAAGAATCTATGGGAATCTCAGTGTGCCTCAAGGGGCCCAATGGCCTGGCATTCTCAGCCTCTT
TGTCAACGTGGCTGGCAGCTCCAACCGTTCTTAGTCGCTCCTCAACCGGGACACTATTACCCGAATC
TCCTACAAAAATGATGCCTACTTTCTTCAAGACCTGAGCCTGGAGCTCCTGTTCCCGGAATCCTTTGGCT
TCATCACCTATCAGGGATCTCTCAGCACTCCTCCCTGCTGAGACTGTACCTGGATCCTCATCGACAG
GGCTCTCAATATTACCTCCCTCCAGATGCACTCCCTGAGACTCCTGAGCCAGAATCCCCCTCCCAGATC
TTCCAGAGCCTCAGTGGTAATGGCCGGCCCTTGACGCCCTGGCCACAGAGCCTTGAGGGGCAACAGGG
ACCCCGGCACCCCGAGAGGCGCTGCCGAGGCCCAACTACCGCCTGCATGTGGATGGTGGTCTCATGG
TCGCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_175708
Insert Size:	987 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_175708.3</u> , <u>NP_783639.1</u>
RefSeq Size:	1346 bp
RefSeq ORF:	987 bp
Locus ID:	308588
Cytogenetics:	1q22