

Product datasheet for RN211720

Tex37 (NM_001109240) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tex37 (NM_001109240) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tex37
Synonyms:	RGD1562720
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN211720 representing NM_001109240 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGGATCGCC

ATGGCACATGTGGTGAACCCCAAAAGTACCAAGTCGACCTAGATATGTACCAAAGCTCATACTGATCG
ACTTCAAGCCCTTCGGGAAGCACAAATACTCCAGAGTCACACCACAAGAGCAGGCCAAGCTTGACGCTCA
GCTCCAAAGCAAAGAGTTTTACCGGCCAAGCCAGCCCAACCCAACTGGAGGATGGGTATCCTGTCT
TTCAAAGAGTCCACATGACTGCCTAGACTTGGGAGTTCTGGCTTCTCCCGCCACGAGACTCAACGG
CTAAGAAGGATGATGGCAGGTTTACCACCACCTGCCATTACGCGTACCCAGCTTCCCTGGCTCTGTACCT
GTCGCAGTACGACCCCTACTGGCTCCACCAGAGGGCAGATTTCCCGTGCCTCATGGAACCAGAGCGGCAG
CCTGCTCCAGAAGTGGGCAAAGGCTACCTTCTACTCCCTGGCTGTCCCTGTGATCATCACCAGAAGGTCA
AGGTCCCATCTTGAATAGGTGGGGCCCTTGATGCCGTTTTACCAGTAG

AGCGGACCGACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-RsrII
ACCN:	NM_001109240
Insert Size:	540 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109240.1, NP_001102710.1</u>
RefSeq Size:	822 bp
RefSeq ORF:	540 bp
Locus ID:	500197
Cytogenetics:	4q31