

Product datasheet for RN211770

Bcl2l10 (NM_053733) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Bcl2l10 (NM_053733) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Bcl2l10
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>RN211770 representing NM_053733 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGTGACCCGCTGCAGGATCGCACTAGACGGCTGCTGACTGACTACATATTGTTCTGCGCACGGGCGC
CGAACACCCTGAGCCACTGCCACGTCTGTTGAGGCGGCCTTGCTGCGCTCTGTGACTAGTCAGATCCA
ACAGGAGCACCAGGATCTTTCAACTCCTCCGCGACTACCAGGGCAACCGCTGGAGCTGGTGACACAG
ATGGCGGATGAGTTGCTCTCCAATGACCAAGAGTTCAACTGGGGCCGCCTGGTGATGCTCCTGGCCTTCG
TGGGGACGCTAATGAACCAAGACAGGACTGTTAAGCGGAGGAGGATCAAAGAAACCGTCTCCTACTGGA
GCGAGACTGCTATCTCATAGTGAGCTTGTGTACAATCGACTCACAGGACGGCATCGCTCCTGGCTGGAG
GCTCACGGTGGCTGGATGGCTTTTGC AATTCTTCAAGAACCCTTACCACCGGCTTCTGGAGAAGAT
TGCTGATCCGGGCTATTCTGTCTGTTTCTTGAACGGCCATCTTTTATATCTGGAAATGTTTA**TAA**

ACGGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_053733
Insert Size:	558 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_053733.1](#), [NP_446185.1](#)

RefSeq Size: 1074 bp

RefSeq ORF: 558 bp

Locus ID: 114552

UniProt ID: [Q99M66](#)

Cytogenetics: 8q24

Gene Summary: Promotes cell survival. Suppresses apoptosis (By similarity).[UniProtKB/Swiss-Prot Function]