

Product datasheet for RN211916

Rhbdl2 (NM_001106684) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rhbdl2 (NM_001106684) Rat Untagged Clone
Tag: Tag Free
Symbol: Rhbdl2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN211916 representing NM_001106684
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGCTGGTCCATGCTGGAGTCCAGCACATCGTGGGAATCTTTTTATGCAGCTCGTCTGGGTATCCCT
 TGGAGATGGTCCACAAAGCCCTCCGAGTGGGCTGGTGTACCTGGCAGGCGTCTTGAGGTTCCCTGGC
 CAGCTCCATTTTGGACCCCTTAAGTCTCTCGTGGTGTCTCAGGAGGAGTCTATGCCTTGATGGGTGGC
 TATTTTATGAATGTTATAGTGAATTTTCGAGAAATGATTCTGCCCTTGAATCGTCAGGCTACTCGTTA
 TCATTCTGATAGTTGCCTCAGACATGGGCTTTGCTCTCTACAGAAGATTCTTTGTCCTGCAAAATGGGTC
 TCCGGTGTCTTCGCTGCTCACATTGCAGGCGGATTTGCTGGGATGTCCGTAGGTTACACTGTGTTACAG
 TGCTTCGACAAGACCCTGCTCAAAGACCCAGGTTCTGGATCTCGATCGCGGCCTATGTAGCTTGTCTCT
 TATTTGCTGTGTTTTTAAACATCTTTCTGTCCCCAGCAAATGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_001106684
Insert Size: 534 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001106684.1](#), [NP_001100154.1](#)

RefSeq Size: 973 bp

RefSeq ORF: 534 bp

Locus ID: 298512

Cytogenetics: 5q36