

Product datasheet for RN211926

Gmppb (NM_001108781) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gmppb (NM_001108781) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Gmppb
Synonyms:	RGD1560458
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN211926 representing NM_001108781 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGGCACTGATCTTGGTGGGTGGGTATGGGACACGGTTGCGGCCGCTGACGCTGAGCACTCCGAAGC
CATTGGTGGACTTCTGCAATAAGCCATTTTGTGCACCAAGTGGAGGCACTTGCCGCGCGGGTGTAGA
CCATGTGATCCTAGCAGTGAGCTACATGTCACAGATGTTGGAGAAGGAAATGAAGGCACAAGAGCAAAGG
CTTGAATCCGAATCTCTATGTCTCATGAAGAGGAGCCTTTGGGAACAGCTGGGCCTCTGGCCTTGCC
GAGACCTACTCTGAGACTGCCGATCCTTTCTCGTTCTCAACAGTGATGTGATATGTGATTTCCCTT
CCAAGCCATGGTGCAATTCATCGACACCATGGCCAAGAGGGCTCCATCTTGGTAACCAAGGTGGAAGAG
CCTTCCAAGTACGGTGTGGTGGTGTGCGAGGCTGACACAGGCCGATTACCGTTTTGTTGAGAAGCCG
AGGATTTGTGTCCAATAAGATCAATGCAGGCATGTACATCCTGAGCCCTGCTGTGCTGCAGCGCATCCA
GCTGAAGCCTACATCCATTGAGAAGGAGATCTTCCAGTTATGGCCAAGGAGGGGCGAGCTATATGCCATG
GAGCTGCAGGGCTTCTGGATGGACATTGGACAGCCTAAGGACTTCTCACTGGCATGTGCCTTCTCTAC
AGTCACTGAGACAGAAGCATCCAGAGAGACTGTACTCAGGCCCTGGCGTTGTGGGAATGTGCTTGTGGA
CCCAAGTGCCCGTATCGGTGAGAAGTGCAGCATTGGCCCAATGTGAGCCTGGTCCCGGTGGTGGTGG
GAGGACGGCGTGTGCATCCGGCGGTGCACGGTGTCTCGTATGCTCACATCCGCTCCCATCTCTGGCTTG
AGTCATGCATCGTGGCTGGCGCTGCCGTGTAGGCCAATGGGTGCGCATGGAGAATGTACAGTGCCTGGG
TGAGGATGTCATAGTTAATGATGAACCTTACCTCAACGGAGCCAGCGTGCTGCCCAACAAGTCTATCGGG
GAGTCGGTGCCAGAACCTCGTATTATCATGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001108781
Insert Size:	1083 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108781.1</u> , <u>NP_001102251.1</u>
RefSeq Size:	1923 bp
RefSeq ORF:	1083 bp
Locus ID:	363145
Cytogenetics:	8q32