

Product datasheet for **RN212022**

Mob3a (NM_001108734) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Mob3a (NM_001108734) Rat Untagged Clone

Tag: Tag Free

Symbol: Mob3a

Synonyms: Mobkl2a; Mobkl2b

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >RN212022 representing NM_001108734
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGTCAAACCCATTCTGAAACAAGTGTTCAACAAGGACAAGACATTAGGCCCAAACGCAAGTTCGAGC
CTGGCAGCGAGCGCTTCGAGCTACACAAGCGCGCGCAGGCATCTCTGAACGCGGGCTGGACCTGCGGT
GGCTGTGCAGCTGCCGCGGGCGAGGATCTCAACGACTGGGTGGCGTGACGTGGTGGACTTCTTCAAC
CGCGTCAACCTCATCTACGTACCATCAGCGATGGCTGTACCGAGCAGTCGTGCCGTGATGTCCGGGG
GCCCAAGTACGAGTACCGCTGGCAGGACGAGCAGCGCTACCGGAAGCCACAGCGCTGTCCGACCCCG
CTACATGGACCTGCTTATGGACTGGATCGAGGTGCAATCAACAACGAGGACCTTTCCCACTAATGTT
GGCACACCGTTTCCCAAGACCTTCTGCAGGCAGTGCGCAAGATCCTGTGCGGGCTCTCCGCGTGTTCG
TGACGTGTACATCCACCATTTGACCGCATCGCACAGATGGGCTCCGAGGCGCACGTGAACACCTGCTA
TAAGCACTTCTACTACTTCGTCACCGAGTTCAGCCTATTGACCCCAAAGAGCTGGAGCCCTGAAGGAG
ATGACATCGCGGATGTGTAC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001108734

Insert Size: 654 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001108734.2</u> , <u>NP_001102204.1</u>
RefSeq Size:	2513 bp
RefSeq ORF:	654 bp
Locus ID:	362833
Cytogenetics:	7q11