

Product datasheet for **RN212076**

Srsf5 (NM_019257) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Srsf5 (NM_019257) Rat Untagged Clone
Tag: Tag Free
Symbol: Srsf5
Synonyms: HRS; Sfrs5; SRp40
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN212076 representing NM_019257
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGTGGCTGTCGAGTGTTCATTGGGAGACTAAATCCCGCAGCGAGGGAGAAAGATGTGGAAAGATTCT
TCAAGGGTTACGGACGGATCAGAGATTTGATTTGAAAAGAGTTTTGGTTTTGTGGAATTTGAGGACCC
AAGGGATGCAGATGACGCTGTTTATGAACCTTGATGGGAAGGAACTTGCAGTGAAAGGGTTACTATTGAA
CATGCTCGGGCTCGGTCTCGAGGTGGAAGAGGTAGAGGACGATACTCCGACCGTTTTAGCAGTCGCAGAC
CTCGAAATGATAGACGAAATGCTCCACCTGTAAGAACAGAAAATCGACTTATAGTTGAGAATTTATCCTC
AAGAGTCAGCTGGCAGGATCTCAAAGATTTTCATGAGACAAGCTGGGGAAGTAACCTTTGCGGATGCACAT
CGACCTAAACTAAATGAAGGGTGTAGTTGAGTTGCCTCTTATGGTGACTTAAAGAATGCTATTGAGAAAC
TTTCTGGAAGGAAATTAATGGGAGAAAATCAAATTAATTGAAGGCAGCAAAAGACACAGGTCAAGAAG
CAGGTCACGATCTCGACAGGAGTTCCTCTAGGTCCCGTAGCCGATCCCGTTACGCAGGAGCAAGTCT
TACAGCCGATCAAGGAGCAGGAGCCGGAGCCGGAGCAAGTCCCGTTCTGGTAGTAGGTCTCCCGTGCCTG
AGAAGAGCCAGAAACGTGGTTCTTCGAGTAGATCTAAATCTCCAGCATCTGTGGATCGCCAGAGGTCCCG
ATCCCGGTCCAGGTCCAGATCAGTTGACAGTGGCAAT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAAATCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_019257
Insert Size: 810 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_019257.3 , NP_062130.2
RefSeq Size:	1501 bp
RefSeq ORF:	810 bp
Locus ID:	29667
UniProt ID:	Q09167
Cytogenetics:	6q24
Gene Summary:	may act as a regulator of alternative pre-mRNA splicing; may play a role in cell cycle regulation [RGD, Feb 2006] Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. All three variants encode the same protein.