

Product datasheet for RN212092

Tas2r139 (NM_001080905) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tas2r139 (NM_001080905) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tas2r139
Synonyms:	Tas2r39
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN212092 representing NM_001080905 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGCAACCCAGCAACTATTGGAAACAAGATTTGCTACCACTGTCCATCTTGATCCTAACACTTGTGG
CCACTGAGTGCACCATAGGTATCATTGCAAGTGGGATCATCACAGTTGTGAATGCAGTGCATGGGTTCA
GAAAAGGGCAGTTCCATAACTACTAGGATTCTGCTTCTCTGAGCGTATCCAGAATAGGCCTCAAAGC
ATCATCTTGATAGAAATGACTTCTCCATATTCAACTTTCTTCTTACAACAGTGTATATATAGAGTCT
CAAGGGTAAGCTTTGTATTCCTAAATTATGTAGCCTCTGGTTTGTCTGCTTTGCTTAGTTTCTTCCACTT
TGTGAAGATTGCAATTTTTCTTACCCCTGTTCTTCAAGCTAAAGTGGAGAATTTCTGAATTGATGCC
TGGCTTCTATGGCTCTCGGTGTTATTTCTTCCAGCTCCAGCATGTTCTTCTGCAATCATAAATACACTG
TGTACAACAACATTTCTCTAAGTAGCAACATCTGCAACTTCACAATGGAAGTCTATGTCGCTGAGGCCAA
TGTGGTCAATGTGGCCTTTTATTAGTTTTGGAATCCTCCCACCTCTGACCATGTTTCAATGCAACAGCT
ACTTCTTAATTTTTCTCTCAGGAGACACCCCTGCACATGAGAAACGGTATGCTGACTCCAGAAATC
CCCGAGTAGAGGCTCATAAGCAGGCCATCAAGGAAACCAGCTGCTTCTCTTCTACATCTTATATGC
AGCTGTTCTGTTCTGTCCACATCCAACATAGCTGATGCCAGTCTCTTCTGGAGTAGTGTCTCAGAATC
AGTCTGCCGGTCTACCCAGCTGGCCACTCAGTTTTACTGATTAGAGCAACCCCTGGCTTAAAAAGAACGT
GGAAGCAACTTCTGTCCAAATCCATCTGCACTTACAAAGTAGATAC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001080905



[View online >](#)

Insert Size:	960 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001080905.1</u> , <u>NP_001074374.1</u>
RefSeq Size:	960 bp
RefSeq ORF:	960 bp
Locus ID:	680188
UniProt ID:	<u>Q67ER9</u>
Cytogenetics:	4q24
Gene Summary:	Putative taste receptor which may play a role in the perception of bitterness. [UniProtKB/Swiss-Prot Function]