

Product datasheet for **RN212109**

Rmdn1 (NM_001031663) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rmdn1 (NM_001031663) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Rmdn1
Synonyms:	Fam82b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN212109 representing NM_001031663 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGCGCTGTCTCGCGACTGTGGCGTCTTCTCCCTCTGCTACGCGGAGCGGAGCAGAGCGTGTGTACG
TGGCGAGGAGCTGGGCGAGCGGTGGGCGCTGTGGCACTTGGGGATTCCGCGGTTCCGAGGTGACAGGAAA
CACAAGAGCCTTCAAGAGAGGCCTTTTGTCTCAGCCTTGTCTTATTTGGGTTTTGAAACCTATCAAGTC
ATTTCTAAGCTGCTGTGGTTCATGCCACAGCCAAAGTTGAAGAAATACTTGAGCAAGCAGACTACCTCT
ATGAAAGTGGAGAAACAGAAAAGCTTTATCAGCTGCTTACACAGTACAAGGAGAGTGAAGATGTAGAGCT
ACTGTGGCGGCTGGCAGCAGCATCACGTGACATAGCACAGCTATGTAAGACTTCAGAAGAAGAAAAGAAA
GTGTTGTTGTATGAAGCTCTAGAGTATGCAAAAAGAGCACTAGAAAAAAGGGGTCAAGCTTTGCTGCTC
ATAAGTGGTATGCAATCTGCATTAGTGATGTTGGGGATTATGAAGGCATCAAGGTTAAGATTGCAATGC
CTATGTCATCAAGGAGCACTTTGAGAAAGCAATTGAACTGAACCCTAAAGATGCTACTTCAATTCACCTT
ATGGGTATTTGGTGCTATACATTTGCTGAAATGCCTTGGTATCAAAGAAGAAATTGCTGAAGTGCTGTTTG
CAATCCTCCTAGTTCACCTATGAGGAGGCCTTAAATACTTTACAGGGCAGAAGAGGTGGATCCAAA
CTTCTACAGCAAAAAGCTGCTGCTTTTAGGAAAGACATACTTGAAGCTAAACAACAAAAGCTTGCTGCA
TTCTGGCTGGTAAAAGCTAAGAGTTACCCAGCACACAGAAGAAGATAAACAGATACAGACAGAAGCTG
CTCAGTTGCTTACAGGTTTG**TGA**

ACGCGTACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001031663



[View online »](#)

Insert Size:	933 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001031663.1</u> , <u>NP_001026833.1</u>
RefSeq Size:	1697 bp
RefSeq ORF:	933 bp
Locus ID:	500419
UniProt ID:	<u>Q4G069</u>
Cytogenetics:	5q13