

Product datasheet for RN212190

Clec5a (NM_001109377) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Clec5a (NM_001109377) Rat Untagged Clone
Tag: Tag Free
Symbol: Clec5a
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN212190 representing NM_001109377
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGAACTGGCACATGACCATCTCAGCTCTTATCGTAGTAGTGATCAAAGTTGTTGGAATGACCTTTTTTC
 TGCTGTACTTCCCACAGTTTTTGGCAAAGTAGTGATGGCTTCATCCCCACGGAGAGCTACGGAACCC
 TAATGTGCAGAATGTCTCACAGATTTTTGGAAGAAGTGATGAAAGCTTCGTTCCCAAGGAGCTTGGA
 ACAGCCTGTCCCAACAACCTGGGATTTTACCAAGGAAGATGTTTTTCTCTACTCTGAATCACCTT
 GGAAAAACAGCAGGGACTATTGTGCAACAAAAGGATCCACACTGGCAATTGTCGACACTCCAAAGAACT
 GAAGTTTCTTCAGGACAGAGCTGGGGCTGAGAATTACTTTATTGGTTTGATACGTCAGCCTGGAGAGAAG
 AAATGGCGCTGGGTCAACAACCTGTGTCAAGGGCAATGTTACCAATCAGGAGCAGAATTCAACTGCG
 TCACTATAGGTCTGACAATGACATTTGATGCTGCATCGTGTGACATCAGCTATCGCTGGATCTGTGAGAC
 CACTCCCAATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001109377
Insert Size: 573 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109377.1, NP_001102847.1</u>
RefSeq Size:	1394 bp
RefSeq ORF:	573 bp
Locus ID:	679787
Cytogenetics:	4q23