

Product datasheet for RN212277

Ltb4r2 (NM_053640) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ltb4r2 (NM_053640) Rat Untagged Clone
Tag: Tag Free
Symbol: Ltb4r2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN212277 representing NM_053640
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGTCTGTCTGCTACCGTCCGCCTGGGAATGAGACGCTGCTGAGTTGGAAGGGCTCGCGGGCCACCGCA
CTGCCTTTCTCCTGCTGGCGGCTTGTGGGACTGCCAGGCAATGGCTTCGTAGTGTGGAGCTTGGCGGG
CTGGCGGCCACCGCGGGCGGCCACTAGCAGCCACACTTGTGCTGCATCTGGCGTAGCCGACGGCGCG
GTGCTGCTGCTCAGCCGCTCTTTGTGGCTTCTGAGCCGACAGGCTTGCCCTGGGCCAGGTGGGCT
GCAAGGCAGTGTACTACGTGTGCGCGCTCAGCATGTACGCCAGCGTGTGCTCACCGGCCTGCTCAGCCT
GCAGCGCTGTCTAGCGGTCACTCGGCCTTTCTGGCTCCCGACTTCGCAGCCGGCCCTGGCCCGCCG
CTGCTGCTGGGGTCTGGCTGGCCGCCCTGGTGTGCGCGTCCCGGCCGGTCTACCGCCACCTCTGGG
GTGATCGCGTGTCAATTGTGCCACCCATCGGCCGTGCACGCTGCAGCTCATCTGAGCCTGGAGACCT
GACTGCCTTCGTCCTGCCTTTTGGGACCGTGTGCGCTGCTACGGCGTGACGCTGGCGGGTTCGCGGGA
GCGCGCTGGGGCTCGGGGCGACAAGGCACGCGGGTGGTGTGCTGCTGGTGGAGCGCCATCGTACTGCCCTTTG
GCTTGCTCTGGGCCCCCTACCACGCGGTCAATCTCCTACAGGCGGTGGCCGCGCTCGCTCCGCCGAAGG
ACCCCTAGCCAGGCTCGGTGGGGCGGGCCAGGCAGCGCGCTGGAACACAGCCTTGGCTTTCTTATG
TCCAGCGTCAACCCGGTCTCTACGTCTTTACTGCGGGTATTGCTGCCGCGGGCGGGCCCTCGGTTCC
TCACTCGACTCTTGAAGGCTCTGGGGAGGCCGAGTAGGCAGCCGCTCTAGGGAGGGTACCATGGAGCT
CCGAACACCCCGAGGCTGAAAGTAGTGGGTCAGGGCAGGGGCTATGGAGACCCTGGAGGTGGGGCAGG
ATGGAGAAAGACAGTCAGGAATGGTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_053640
Insert Size: 1077 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_053640.1](#), [NP_446092.1](#)

RefSeq Size: 1077 bp

RefSeq ORF: 1077 bp

Locus ID: 114098

UniProt ID: [Q924U0](#)

Cytogenetics: 15p13

Gene Summary: mouse homolog is a chemoattractant for myeloid leukocytes [RGD, Feb 2006]