

Product datasheet for **RN212417**

Zbtb33 (NM_001109314) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zbtb33 (NM_001109314) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Zbtb33
Synonyms:	RGD1566309
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001109314.1](#), [NP_001102784.1](#)

RefSeq Size: 2561 bp

RefSeq ORF: 2013 bp

Locus ID: 501506

Cytogenetics: Xq35