

Product datasheet for **RN212440**

Cpm (NM_001108098) Rat Untagged Clone

Product data:

| | |
|----------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Cpm (NM_001108098) Rat Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Cpm |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Cell Selection: | Neomycin |
| Fully Sequenced ORF: | >RN212440 representing NM_001108098 Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACGGCGCGCCTGTGGCTGGGGCTTCTGCTGCCCGTGGTGGCCGCTTGGATTTCCGATACCACC
ACCAGGAAGGGATGGAGGCTTTCCTAAAGAGCGTGGCTCAGAACTACAGTTCGATCACGCACTTGCACAG
CATCGGAAATCGGTGAGAGGTAGAACTGTGGTCTTGTGTGGGGCAGTCTCAAAGGAGCACAGA
GTCCGGATTCCAGAGTTCAAATACGTGGCTAATATGCACGGCGATGAGACTGTGGGCGAGAAGCTGCTGC
TCCACTTGATCGACTATCTGGTAACCAAGTTATGGGAAGGACGCCGAAATGACACGCCTGATCGACAGCAC
CCGGATCCACATCATGCCTTCCATGAACCCAGATGGATTCGAAGCCGTCAAGAAGCCGACTGTTACTAC
AGTAATGGAAGGGAAAATAACAATAATTATGACCTGAACAGAACTTCCCTGATGCCTTTGAAAATAACA
GTGTAACCAAGCAGCCGAGACCCTGGCGGTATGACGTGGTGCAGAACCCGAGACATTTGCCTCTCTGC
GAATCTCCATGGGGTGCCCTGGTGGCCAGTTACCCCTTGATAATGGCGTGCAAGCCACGGGGACTG
TTGTCCCGAAGCCTGACTCCAGATGACGATGTTTTTCAACACCTGGCGTACACCTATGCCTCACGGAATC
CCAACATGACAAAAGGGGATCAGTGTAAAAACAAGAGGAGCTTCCCAACGGGATCATCAACGGGACTC
CTGGTATCCACTGCAAGGGGAAATGCAAGATTACAACATACATCTGGGCCAGTGCTTTGAAATCACCTG
GAGCTGTCTGTGTAATAACCTCGGGAGGAGAAGCTACCGTCTTTTGGAAATGATAACAAAGCCTCTT
TGATCGAGTACATAAAACAGGTGCACCTAGGTGTAAAGGGTCAAGTGTGTTGATCAGAGTGGAAATCCATT
GCCGAATGTAATTGTGGAAGCCCAAGACAGAAAACATACTGCCCGTTCAGAACCAACAGCTTGGAGAG
TATTACCTGCTTCTTTGCTGGGTCTACGTGATAAATGTTACAGTCCCTGGACATGACTCACACCTCA
CGAAGCTTACAATTCCATCGAAATCCAGCCCTCAGTGTCTTAAAAAGGATTTTACCTCCCGTTGCG
ATGGCAGCCAGATTCGTCTCCGTATCAAATCCTTCGTGCCAATGATTCCTCTGTACAAATTAAGCCA
AGCCACTCTGTGCAACAAAGCCTAGTTGTTTGTATTTTCGTGACTCTTTTGCACATATTTTCAAAT
AA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

| | |
|-------------------------------|---|
| Restriction Sites: | Sgfl-Mlul |
| ACCN: | NM_001108098 |
| Insert Size: | 1332 bp |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001108098.1, NP_001101568.1</u> |
| RefSeq Size: | 5450 bp |
| RefSeq ORF: | 1332 bp |
| Locus ID: | 314855 |
| Cytogenetics: | 7q22 |