

## Product datasheet for **RN212467**

### **Rpa2 (NM\_021582) Rat Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rpa2 (NM\_021582) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Rpa2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN212467 representing NM\_021582  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**

ATGTGGAATAGCGGATTCGAAAGCTTCAGCAGCTCCTCCTATGGCGCAGCGGGCGGCTACACGCAGTCCC  
CCGGTGGCTTCGGGTCACCCACACCGTCGCAGGCAGAGAAGAAATCAAGAGCCCAGCCCAGCACATTGT  
GCCCTGTACCATATCTCAGCTGCTCTCTGCTACTCTGACTGATGAGGTGTTCAAATTTGGGGATGTCGAG  
ATTTACAGGTCACTATTGTGGGAATAATCAGACATGCAGAGAAAGCGCCAACCAACATTGTCTACAAAA  
TAGACGACATGACCGCAGCGCCCATGGACGTTCCGACAGTGGGTTGACACAGATGATACCAGTGGCGAGAA  
CACTGTGGTTCCTCCAGAAACATACGTGAAAGTTGCTGGCCACCTCAGGTCTTTTCAGAACAAAAAGAGC  
CTGGTGGCCTTTAAGATCATTCTCTGGAGGATATGAATGAGTTCACGGCGCACATTCTGGAAGTGGTCA  
ATTCACACTTGATGCTCAGCAAAGCCAACAGCCAGGCCTCTGTAGGGAGACCATCTATGAGCAACCCAGG  
AATGGGTGAGCCAGGAACTTCAGTGGGAATAACTTCATGCCAGCAAATGGCCTCACCGTGGTCCAAAAC  
CAGGTGCTGAATTTGATTAAGGCTTGCCCAAGACCTGAAGGATTGAACTCCAGGATCTCAGGAGCCAGC  
TCCAGCACATGCCTGTAGCCTCAATCAAACAAGCAGTGGATTTTCTATGCAACGAGGGCCACATTTACTC  
TACAGTGGATGATGATCACTTTAAATCCACGGACGCAGAT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_021582  
**Insert Size:** 813 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_021582.1</a> , <a href="#">NP_067593.1</a>
<b>RefSeq Size:</b>	1436 bp
<b>RefSeq ORF:</b>	813 bp
<b>Locus ID:</b>	59102
<b>UniProt ID:</b>	<a href="#">Q63528</a>
<b>Cytogenetics:</b>	5q36
<b>Gene Summary:</b>	may bind uracil-DNA glycosylase; may play a role in base-excision repair [RGD, Feb 2006]