

Product datasheet for RN212497

Pdik1l (NM_001107984) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pdik1l (NM_001107984) Rat Untagged Clone
Tag: Tag Free
Symbol: Pdik1l
Synonyms: RGD1307476; Stk35l2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN212497 representing NM_001107984
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGTGAAGTACGCAAAAGTACGATCTAATACGGGAGGTAGGCCGAGGTAGTTACGGTGTGTGTATG
AAGCAGTCATCAGAAAGACCTCTGCGAGAGTAGCCGTGAAGAAAATTCGATGCCACGCACCCGAAAATGT
TGAAGTACCCCTTCGAGAGTTCTGGGCACTAAGCAGTATCAAGAGCCAACACCCAATGTGATTCACCTTG
GAAGAGTGCATCCTACAGAAAGATGGGATGGTGCAAAAGATGTCCACGGCTCTAATTCTCCCTTTATT
TACAGCTTGTAGAGACTTCATTAAGGAGAGATCGCCTTTGACCCAGGAGCGCCTACTATTTGTGGTT
TGTGATGGATTTTGTGATGGAGGAGATATGAATGAGTATCTGTTGTCCAGGAAACCCAATCGAAAGACT
AACACCAGCTTCATGCTTCAGCTCAGCAGTGTCTGGCCTTTTGCACAAAAACCAAGATAATCCACCGAG
ACCTGAAGCCTGATAACATCCTGATTTCTCAGAGCAGATTGGATACCAGTGACTTGGAAACCCACTGAA
AGTGGCTGATTTTGGCCTAAGTAAAGTTTGCTCAGCATCTGGGCAGAACCAGAAAGAACCTGTCAGTGTA
ACAAGTGTTCCTTTCCACGGCATGTGGAACGGACTTCTACATGGCTCCGGAAGTATGGGAAGGACACT
ACACGGCGAAAGCTGACATTTTGTCTGGGAATTATCATCTGGGCCATGCTGGAGAGGATTACCTTCAT
AGACACAGAGACGAAGAAGGAGCTCTGGGGAGTTACGTAAGCAGGGGACCGAAATCGTGCCTGTCGGG
GAGGCACTTCTGAAAAACCCAAAATGGAACCTTCTCATCCTGTGAAGAAAAAGTCTATGAATGGCGAA
TGAAACAAGTATAAGGAAATGCTGGCCGAAACCCCTCAGGATCGCCGGATGCTTTTGAAGTAACT
CAGATTAGTCCAAATTGCATTTAAAGACAGCAGCTGGGAAACGTGA

ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAAGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001107984



[View online »](#)

| | |
|-------------------------------|---|
| Insert Size: | 1026 bp |
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001107984.1</u> , <u>NP_001101454.1</u> |
| RefSeq Size: | 4489 bp |
| RefSeq ORF: | 1026 bp |
| Locus ID: | 313609 |
| Cytogenetics: | 5q36 |