

Product datasheet for **RN212501**

Tial1 (NM_001013193) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tial1 (NM_001013193) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Tial1
Synonyms:	MGC109552
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN212501 representing NM_001013193 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGATGGAAGACGACGGACAGCCCCGGACTCTATACGTAGGTAACCTCTCCAGAGATGTGACAGAAGTCC
TTATCCTTCAGTTATTCAGTCAGATCGGACCCTGTAAAAGCTGTAAAATGATAACAGAGCAACCCGATAG
CAGAAGGGTCAACTCTTCTGTTGGATTTCTGTTTTGCAGCATAACAAGCAATGACCCATATTGCTTTGTG
GAATTTTATGAACACAGAGATGCAGCTGCTGCATTAGCTGCTATGAATGGGAGAAAAATTTGGGAAAGG
AGGTCAAAGTAACTGGGCAACAACACCAAGTAGCCAGAAAAAGACACTCCAATCACTCCATGTGTT
TGTTGGGGATTTGAGTCCAGAAATTACAACAGAAGATATCAAGTCAGCATTGGCCCTTTGGTAAATA
TCGGATGCCCGGGTAGTTAAAGATATGGCAACTGGAAAGTCCAAAGGCTATGGTTTTGTATCTTTTTATA
ACAAACTGGATGCAGAGAATGCAATTGTGCATATGGGGGGTCAAGTGGTGGGAGGTCGACAAATCAGAAC
CAATTGGGCCACACGTAAACCACCTGCCCTAAAAGTACACAAGAACTAACACTAAGCAGTTGAGATTT
GAAGATGTAGTAAACCAGTCAAGTCCAAAAACTGCACCTGTGTACTGTGGGGGGATCGCTTCGGGGCTGA
CAGATCAACTTATGAGACAAACATTTTCACCATTTGGACAATATGGAGATAAGAGTTTTTCCAGAGAA
GGGGTATTCGTTTGTGAGATTTTCAACTCATGAAAGTGCAGCCCATGCTATTGTTTCTGTGAACGGTACT
ACAATCGAAGGACATGTGGTAAATGCTATTGGGGTAAAGAATCTCCTGATATGACTAAAAACTTCCAAC
AGGTTGACTACAGTCAGTGGGGCCAATGGAGCCAAGTGTATGAAACCCACAACAGTACGGACAGTATAT
GGCAAAATGGGTGGCAAGTACCTCCCTATGGAGTATATGGGCAGCCATGGAACCAACAAGGATTTGGAGTA
GATCAATCACCTTCTGCTGCTTGGATGGGTGGATTTGGTGCTCAGCCTCCACAAGGACAAGCTCCGCCCC
CTGTCATACCTCTCAAACCAAGCTGGGTATGGAATGGCAAGTTTCCAACACAGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001013193
Insert Size:	1179 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001013193.1</u> , <u>NP_001013211.1</u>
RefSeq Size:	1680 bp
RefSeq ORF:	1179 bp
Locus ID:	361655
Cytogenetics:	1q37