

Product datasheet for **RN212524**

Fbxl13 (NM_001103355) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fbxl13 (NM_001103355) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Fbxl13
Synonyms:	MGC95282
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF:	>RN212524 representing NM_001103355 Red=Cloning site Blue=ORF Orange=Stop codon TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCCCGGATCGCC ATGGCTTCACTACGCAACGCCCCCCCCAGGCTAAGGAACACTTCAAGGACAACACATTCCTCAAATCT GTGAGGCGCTGTTATGTGGTTTACTTGTGACATGCCCTGAGGACCCATTGAAATATTTAGAGAGCATGAT CCTAGCTATCATTGAAAATGGTCTGAAAACCTTCTTTGGGACATGTGCATTACCCCATCACTGAAGCCA AAAGTTGGAAGACTGTCTGAACTTACTTAGATGAGCTTTTGGACTGGATGACCAGCTGATGACTCCAG AGTTGATGATAAAGCATGCACATTTTACACTGGCCATTTAGTAAAGAGCCACTTCGGCGTTTGAAGGA AATTGCAATACCTCATGCAAATCAAGAGGCTATCATGGCTGAGAAAATGGACAAGGCAATAGCCACGAC AATTTGAGATGCCAAAAACACATTTTAAATCGTGGTTCTCTACACAGTTAAGAGAAGAGAAAGGCTTA TAGCTACACTGCTCCGACTGCGCCATTTGTTCTATACGCAAAAAACAAAGATTATCCTAGCAAAATGGAA GGAGAAAGCCAGGCACAAGAGCAAAAAAAGAGAAGATGACTTGATTACAAAACATGAATTCAGTTGAAA AAATGGAAGTTCAAGCTATCCTTGGAGAAACCGATTAAACATTGAGGGAAGCTTCTCCGATTTCCGGGTGA CGGACAAGCGCATAGTGTGTTGACATTTCACTGCTACCCGACCAAGCAGTAGTCCAGATTTTCGTGTACCT CACCTTTAAAGACTTGGTCTCGTGACGTCAGTCAACCGCTCCTGGATGTCGATGATTACAGAGAGCTCC TTGTGGAATTCATTGATTTTCCACAGTAAAAACATTCAGATAGATGTGTGGTGACCACTTTACAAA AGTGGCGTCTAAATGTCTGCGTTTAACTTCCGAGGCTGCGTCTTACGAGCGAAAACCTGAAATCTGT CAGCCACTGTAAAACTTGCAAGAGTTGAATGTCTCCGACTGCCGCTTTTACAGATGAATCCATGAGA CACATTTCCGAGGGCTGTCTGGGGTCTGTATCTCAATCTATCCAACACAACCATCACCAACAGGACAA TGGCAGTCTGCCAGGTAATTCACAACCTCCAGAACCTGAGTCTGGCTTACTGCCGAAAGTTACCCGA CAAAGGCTTACAATACCTGAACCTGGGGAATGGATGCCACAAGCTCATCTATCTGGACCTTCCGGCTGC ACCCAGATCTCAGTCCAAGGCTTCCGGAACATTGCGAACAGCTGCACGGGAATCATGCACCTGACCATCA ATGACATGCCAACGCTGACGGACAACCTGTGTGAAAGTTTTGGTTGAGAAGTGCTCTGATTTTCATCGGT GGTTTTTATTGGCTCCCAACACATCTCGGACTGTGCTTCAAAGCTCTTCTGCCTGTGACCTCAAAAAG ATCCGATTGGAAGGAAACAAAAGGATCACCGACGCGTGCTTAAAGTCGGTAGACAGGAACTACCCAGGCA TCAGTCACATTTACATGGTGGACTGCAAGGGGCTGACAGATAGCAGCCTCAAGTCACTTTCAGTTTTGAA GCAATTGACGGTGTGAACTTGACAAATTGCGTCAGGATTGGCGATATCGGCCTGAGGCAATTTTTCGAC GGTCTGCCAGTGTGAACTAAGAGAGCTCAACTTGGTAAGTCTCCCTGCTGGGCGACACCTCTGTCA TAAGACTGTGAGAGCGATGCCCTAATTTACACTACTTGAATTTACGAAATTTGTGAACATCTGACTGACTT GGCAATTGAATATATTGAAGCATGTTATCCTTGATATCAATAGATCTTCTGGAACACTCATCTCCAAT GAGGCAAGCCATTGA AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC TGGATTACAAGGATGACGACGATAAGGTTTAA
Restriction Sites:	Sgfl-RsrII
ACCN:	NM_001103355
Insert Size:	1905 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001103355.1, NP_001096825.2

RefSeq Size: 2256 bp

RefSeq ORF: 1905 bp

Locus ID: 100125363

Cytogenetics: 4q11