

Product datasheet for **RN212543**

PPP1CB (NM_013065) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: PPP1CB (NM_013065) Rat Untagged Clone
Tag: Tag Free
Symbol: PPP1CB
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN212543 representing NM_013065
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGCGGACGGGGAGCTGAACGTGGACAGCCTCATCACCCGCTGCTGGAGGTACGAGGATGTCGTCGG
GAAAAATTGTGCAGATGACTGAAGCAGAAGTCCGAGGACTGTGTATCAAGTCTCGTAAAATCTTTCTTAG
CCAGCCTATTCTTTGGAATTGGAAGCGCCACTGAAGATTTGTGGAGATATTCATGGACAGTATACAGAC
TTACTGAGATTATTTGAATATGGAGTTTTCCACCAGAAGCCAATATCTTTTCTTAGGAGATTATGTGG
ACAGAGGAAAGCAGTCTTTGAAACCATCTGTTTGTATTGGCTTACAAAATCAAATACCCAGAGAAGT
CTTTCTTCTACGAGGAAACCATGAGTGTGCTAGCATCAACCGCATTTATGGATTCTATGATGAGTGCAAA
CGAAGATTTAATATTAATTGTGGAAGACATTCAGTATTGTTTAAATTGTCTGCCTATAGCTGCTATTG
TTGATGAGAAAATCTTCTGCTGTCATGGAGGACTGTCACCAGACCTACAGTCTATGGAACAGATTCGGAG
AATTATGAGACCCACTGACGTACCTGATACAGGTTTGCTTTGTGATTTACTGTGGTCCGACCCAGATAAG
GATGTACAAGGCTGGGGAGAAAATGATCGTGGTGTTCCTTTACTTTGGAGCTGATGTAGTCAGTAAAT
TTCTGAATCGTCATGATTTGGACTTGATTTGTCGAGCCCATCAGGTGGTGAAGATGGATATGAATTTTT
TGCTAAACGACAATTGGTAACTTTATTTCTGCCCAAATTAAGTGGCGGAGTTTGACAATGCTGGTGGT
ATGATGAGTGTGGATGAACTTTGATGTGTTTCATTCCAGATATTGAAACCATCTGAAAAGAAAGCTAAGT
ACCAATGATGGTGGCTGAATTTCTGGACGTCTGTCACTCCGCCTCGAACAGCTAATCCACCGAAGAAAAG
GTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_013065
Insert Size: 984 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_013065.2 , NP_037197.1
RefSeq Size:	2942 bp
RefSeq ORF:	984 bp
Locus ID:	25594
UniProt ID:	P62142
Cytogenetics:	6q13
Gene Summary:	one of several neuronal isoforms of protein phosphatase 1 [RGD, Feb 2006]