

## Product datasheet for RN212557

### Ggps1 (NM\_001007626) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ggps1 (NM\_001007626) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Ggps1  
**Synonyms:** Crlf3; rGGPS1a; rGGPS1a1; rGGPS1a2; rGGPS1a3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN212557 representing NM\_001007626  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGAGAAAACAAAGAGAAAGCTGAAAGAATCTTCTAGAGCCTTATAAGTACTTACTTCAATTACCAG  
GTAACAGGTGAGAACCAAACTTTCACAGGCATTTAATCACTGGCTGAAAGTTCCAGAAGACAAGCTACA  
GATTATCATTGAAGTACTGAAATGTTGCATAATGCCAGTTTACTCATTGATGATATTGAAGACAGTTCA  
AAGCTCCGACGTGGTTTTCCAGTTGCTCACAGCATCTATGGTGTCCCATCTGTCATTAACCTGCCAATT  
ATGCTACTTCTTGGCCTGGAAAAAGTCTTAACACTTGATCACCCAGATGCAGTGAAGCTTTTCACACG  
ACAGCTTTTGAACCTTATCAGGGACAAGGCCTAGATATTTACTGGAGAGATACCTACACTTGTCCGACT  
GAAGAAGAATAAAGCCATGGTGTGCAGAAGACAGGTGGCTTGTGGATTAGCCGTAGGCCTCATGC  
AGTTGTTCTCTGACTACAAAGAAGATCTAAAGCCACTGCTGGACACACTTGGGCTCTTCTCCAGATTAG  
AGATGATTACGCCAATCTACACTCAAAGAATACAGTAAAACAAAAGTTTCTGTGAAGACCTGACAGAA  
GGAAAGTTCTCATTTCCTACTATTCATGCCATTTGGTCAAGGCCAGAAAGCACCCAGGTACAGAATCC  
TGCGCCAAAGAACAGAAAATATAGACATTAAGAAAGTACTGTGTGACAGTACCTGGAGGATGTAGGTTCTTT  
TGAATATACGCGATACACTCTTAGAGAGCTCGAAGCTAAAGCCTACAAACAAATCGAGGCCTGCGGTGGG  
AACCTTCACTAGTGGCTTTAGTCAAGCATTTAAGTAAGATGTTACAGAAGAAAATGAAT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001007626  
**Insert Size:** 903 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001007626.2</a> , <a href="#">NP_001007627.1</a>
<b>RefSeq Size:</b>	2620 bp
<b>RefSeq ORF:</b>	903 bp
<b>Locus ID:</b>	291211
<b>UniProt ID:</b>	<a href="#">Q6F596</a>
<b>Cytogenetics:</b>	17q12.1
<b>Gene Summary:</b>	Catalyzes the trans-addition of the three molecules of IPP onto DMAPP to form geranylgeranyl pyrophosphate, an important precursor of carotenoids and geranylated proteins.[UniProtKB/Swiss-Prot Function]