

## Product datasheet for **RN212655**

### RT1-T24-4 (NM\_001008826) Rat Untagged Clone

#### Product data:

|                           |                                                                                    |
|---------------------------|------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                |
| Product Name:             | RT1-T24-4 (NM_001008826) Rat Untagged Clone                                        |
| Tag:                      | Tag Free                                                                           |
| Symbol:                   | RT1-T24-4                                                                          |
| Synonyms:                 | H2-T24; RT1-149                                                                    |
| Mammalian Cell Selection: | Neomycin                                                                           |
| Vector:                   | pCMV6-Entry (PS100001)                                                             |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                               |
| Fully Sequenced ORF:      | >RN212655 representing NM_001008826<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGTGGGCTCTGATCTTCTGCTTCCTGAGTCACCCCCGAGACACTGGAGCCCGATCACACTCTCTGCACT  
ATTCTTACTCAGCCGTGACAGAGCCGGGCCAGGAGTTCTGCATTTTGGCAGTGCGTTGTAGATAA  
CCAGCCCTTCATCCACTATGACAGTAGGAAGATGGAGGCAGAGCCTTGCTGCTGACTGGCTGAGGGAAAAT  
CCATATTACTTTGATGAGGAGACCCAGGTTTTACCAATCGGATGAAGATCTTCCAGTTGAGCCTGAGAA  
ACATTGAGCAATACTACAACAGCTCTGGTGACCGGAGTCAGAGAGCAGATGGCTTCCACCAGCAAGCAGG  
CCCTCACACCCTCCAGTTCACGTACGGCTGTGAGATGCAGGACAACAGGACCACAGGGCACTGGCAGTAT  
GGCTACGATGGCAGGGATTACCTGACCTTGGACTTGGACTCTATGCAATATACAGCAGCCACCTTCATTG  
CTGGCTACACCAAGCACAAAGTGGCAAAACAATGAGTACTGGCTGGAGAGAGAGAAGACCTACTTGGAGAA  
GGAGTGTGTCTGTGGCTAAAGAGATACTTAACATATGGGAGGAGAGAAGTTTACCCGAAGTACCCCTCCA  
CAGACAACAGTGACACACCACCCAGACCTGAAGGTGATGTACCCCTGAGGTGCTGGGCCCTGGGCTTCT  
ACCCTGCTGACATCACCTGACCTGGCAGTTGAATGGGGAGGACCTGACCCAGGAAATGGAGCTTGTGGA  
GACCAGGCCTGCAGGGGATGGAACCTTCCAGAAGTGGGCATCTGTGGTGGTGCCTCTTGGGAAGGAGCAG  
AAGTACACATGCCATGTGGAGCATGAGGGGCTGCCTGAGCCGCTTACCCAGAGATGGGAGCCTCCTTCAT  
CCACTGATTCCAACATGGAACCAATGTCAATTTATGTTGTGCTTGGAGCCGCTCTGCTGTGGCCATTAT  
TGGAGCTGTGATCATTGGAGCTGTGGTGGCTGTTGTGAGGAAGAGGAGGAGAAACACAGGTGGGAAAGGA  
GTCTACGCCCTGCTCCAGGCAGGGACAGCTGCCAGAGCTCTGATGTGTCTCTCCAGATTGTGAAGCAT  
GA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Restriction Sites:</b>     | Sgfl-Mlul                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>ACCN:</b>                  | NM_001008826                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Insert Size:</b>           | 1122 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                          |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u>NM_001008826.3, NP_001008826.2</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>RefSeq Size:</b>           | 1678 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 1122 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Locus ID:</b>              | 414784                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 20p12                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |