

## Product datasheet for **RN212725**

### **Gpr150 (NM\_001109173) Rat Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Gpr150 (NM_001109173) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Gpr150
Synonyms:	RGD1564729
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >RN212725 representing NM\_001109173  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGAAGATCCCTTCAGCCCTCAACTCTGTCTCCGGTGTCTAACCTATCCCTACCCATTGGCAGAGCT  
 GGGGTCTCAACCTCACTTCTAAGCAGGGAACCTCTGTGCCGGGACCGCAGCCGCCTCCACGCTGGCCACC  
 TAGCCACCGGATCCATCTGGTCTTCTTGGGGATCATCCTGGCGGCCTCGGTAGCTGGCAACACCACAGTG  
 CTCTGCCGCCTGTGCGGCAGTAGCAGCCGACCCTGGCTGGGTCCCAAGCGTCGCAAGATGGACTTCTGCTG  
 TGGTGCAGCTAGCAGCGGCTGACCTGTACGCGAGTGGAGGCACCGCGCTGTGCGAGCTGGCTTGGGAGTT  
 GCTGGGCGACCCGCGACCTGCTCTGGGCAACTGGCTTGGCGCTTCCAGTTGCTGCAGACATCTGGG  
 CGGGGCGCCTCCGCCACCTGGTGGCGCTATTGCTTTGGAGCGCAGCTCGCAGTGCAGGTCCTCAGG  
 GCCACAGTTGCCCGCACGCGCCCTGACTACCCTGAGCTGGCTGTGGCTCTGCTGCTTGTCTGCCGCC  
 CACCTTCTGGTGGTGGGAGCGCCTCCATCTCAACCGCCAGCGCCTGGCCCGGGAGCATCGCTGC  
 CAGGACATCTTCGCCCCCTGCCGCGCTGGCACCTGCAGGTCTACGCGCTCTACGAGGCTATCGTGGGCT  
 TCGCGGCGCCAGTCGCGATCTTGGGTGTCTTTCGCGTCACTGTGTGTCTGGTGGCGGCGGGGTC  
 TCAGGCTTCGGTGGTGGGATGCCCTGGTCCGCCAGTATGGCTCGAGCCTCCTTACCCAGTGCAGTGGCC  
 CGCGCCAAGGTGCAGAGCTTAAAATGAGCCTGGCACTGGCACTGCTCTTCGTGGGCTGCAGCTGCCCT  
 ACTTTGCAGCTCGCTTGGCAGCCGCTGGTCTCGGGGACCGGGGTGATTGGGAGGAGAGAGCCTGGT  
 AGCCGCGATGCGGGTGTAGAGGTGGCCAACAGCGCTGTCAATCCTTTGGTCTACCTTTTCTCCAGACG  
 GGTGACTGCCGTCTCTGGCGCGGCTACGGAGGCGTATAGGAGTCTGTGCTGCGTGCAGGAGGAAGAAG  
 CCGAGATCAGCGGAGGGGCTGGGACCAACAGGCACTCCACCGCCATCGTGGCCCCACCCCAATTATCA  
 CCATGCCCGGCGGAGGAGCGGGACCAAGCGGGCTGCGCCACCCCGCCAGCCCCAGACCTCTGCC  
 TGCTCCTGTGAAAGTGCCTT**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja1898\\_b10.zip](https://cdn.origene.com/chromatograms/ja1898_b10.zip)

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001109173

**Insert Size:** 1284 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001109173.1, NP_001102643.1</u>
<b>RefSeq Size:</b>	2137 bp
<b>RefSeq ORF:</b>	1284 bp
<b>Locus ID:</b>	499486
<b>Cytogenetics:</b>	2q11