

Product datasheet for **RN212754**

Csf3 (NM_017104) Rat Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Csf3 (NM_017104) Rat Untagged Clone
Tag: Tag Free
Symbol: Csf3
Synonyms: Gcsf
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >RN212754 representing NM_017104
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAAGCTGATGGCCCTACAGCTGCTGCTGTGGCACAGTGCCTATGGTCAGGACAAGAGGCCATCCCC
TGCTCACCGTCAGCTCTCTGCCACCATCCCTGCCTTTGCCCGAAGCTTTCTGCTTAAGTCTTGGAGCA
AGTGAGGAAGATTCAGGCCAGGAACACTGAACTGCTGGAGCAGCTGTGTGCCACCTACAAGCTGTGTCAT
CCGGAGGAGCTGGTGCTGTTGGCCACTCTCTGGGTATCCCAAGGCTTCCCTGAGCAGCTGCTCCAGCC
AGGCCCTTGACAGACAAAGTGCCTGAGCAACTGCACAGCGGGCTTCTCCTCTACCAAGGTCTCCTGCA
GGCCCTAGCAGGCATTTCTCTGAGCTGGCCCCACCTTGGATATGCTTCACCTGGATGTCGACAACTTT
GCCACCACCATCTGGCAGCAGATGGAAAGCCTCGGGGTGGCCCCACTGTGCAGCCCACCCAGAGCACCA
TGCCAATCTTCACTTCTGCCTTCCAGCGTCGGGCAGGAGTCTGGTCACTTCATACCTGCAGAGCTT
CCTGGAGACAGCTCACCACGCTCTGCACCACTTGGCC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_017104
Insert Size: 600 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_017104.2</u> , <u>NP_058800.2</u>
RefSeq Size:	1297 bp
RefSeq ORF:	600 bp
Locus ID:	25610
Cytogenetics:	10q31
Gene Summary:	putative hematopoietic growth factor for neutrophils; used as a treatment for cyclic hematopoiesis in humans and dogs [RGD, Feb 2006]