

## Product datasheet for RN212781

### Tmie (NM\_001109299) Rat Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tmie (NM\_001109299) Rat Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Tmie  
**Synonyms:** RGD1562523  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >RN212781 representing NM\_001109299  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGGCCGGGAGGCAGCGTGGCTCGGGCGGCTCTGGGCGCTGGGCGGCCGCGCTGGGGCTTGCCTGG  
 CGGGTGTGCGCTCGCAGCTGGTAGAGCCAGCACGGCCCCGCCAAGCCCAAGCCGCCCCACTGACCAA  
 GGAGACCGTGGTGTCTGGGACATGCGCCTGTGGCATGTGGTGGGCATCTTTCCCTCTCGTGCTGCC  
 ATCATCATCAGCTGTGCTGTGTGTTCAACTGCCGTGTACCGCGGACACGGAAGGAGATCGAGGCTCGGT  
 ATCTACAGCGAAAAGCAGCCAAGATGTACACAGACAAGCTGGAGACCGTGCCCTCCCTCAACGAACCTCAC  
 AGAAATCCCTGGAGAGGACAAGAAGAAAAGAAGAAGGACAGTGTGGACACCGTGCTATCAAGGTAGAG  
 GAGGACGAGAAGAATGAGGCAAAGAAGAAGGAGAGAAGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_001109299  
**Insert Size:** 462 bp  
**OTI Disclaimer:**

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001109299.1, NP_001102769.1</u>
<b>RefSeq Size:</b>	2599 bp
<b>RefSeq ORF:</b>	462 bp
<b>Locus ID:</b>	501061
<b>Cytogenetics:</b>	8q32