

Product datasheet for **RN212790**

Vom1r22 (NM_001008900) Rat Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Vom1r22 (NM_001008900) Rat Untagged Clone
Tag:	Tag Free
Symbol:	Vom1r22
Synonyms:	V1rf2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>RN212790 representing NM_001008900 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTATGTATGTATGATTTGAAGGCTGGCCTCATCTTCCCTGACAGTGTGGCTGAGAGCTGGATGCA
 ACATCATGGCAGCCAGTGAAGTGGCTGTAGGAGGGATCTTCTGTTACAGACTGTGGTTGGGTTTCTGGG
 CAATTGTTTTCTTCCACCATTATCTATTGCTTAACCTCAATGGGTTTAGGCTAAGGTCCAGATACTGG
 ATTGTGATACACTTGAGTATAGCCAACATCTTAATTCTGTTGTGTAAAGGAGTGCCACAGACAATGTTGT
 CTTTTGGGTTAAAGACTTCCTCAGTGACTTTGGATGCAAAGGGACTTACTATCTCCACAGAGTTGGCAG
 GGGTGTGTCCATAGGCAGCACTACCTTCTGAGTGTCTACCAAGCCATCACTATTAGCCATTGGAATCC
 AGATGTTCAAGTTTAAAGGAAGACGACAAGAACAATTGGTTCTTCTGTTTACCTAAGTTGGATTCTGT
 ACATGCTTATCAGTAGCGTTAACCTTGTGTATGTGAGAGCAAAATACATCAAGAACGACACAGCAAAAT
 AATAGATATGGGATACTGTGTTACTACTGCAGATAAAATCAGAGACTTACTGTATGTAGCACTGTTTCA
 GTTCTGATTTTCTTTTATGGGTTTATGGTCTGGGCTAGCACCATTACGGTTCTCATACTTCATAGGC
 ACAAGCAGAAAATGCAACACATGTCCAGTGCCAAAGTATCTTCCAGATCATCTCCTGAGACAAAAGCTGT
 CCAAACCATCCTTCTCCTGGTGAGCACTTTGTCTTTTACACACTTTCCTGCATCTTTATTGTATGT
 ATAACCTTTTCTAGTCATCACAATGGCTGTTGGTGAATGTTTTGTATTACTTTCTGGATGTTTCCCAA
 CTGTCAGCCCCTTCTTAATCTTAAGACCTGAGTCTACTGTGTTTATGAGTTTTCTTTTATCTCTACAAA
 AAAAAAGAAACAAGCAGGGGATATTATCT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001008900
Insert Size:	1014 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001008900.1, NP_001008900.1</u>
RefSeq Size:	1014 bp
RefSeq ORF:	1014 bp
Locus ID:	494252
Cytogenetics:	1q12